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# Laryngeal Tube LTS-D

---

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
## INTENDED USE

The LTS-D is a supraglottic airway with drain tube.

## INDICATIONS

- As alternative adjunct to temporarily secure the airway during difficult airway management and cardiopulmonary resuscitation in clinical as well as preclinical environments (current guidelines of the European Resuscitation Council, American Heart Association and American Society of Anesthesiologists).
- For mechanical or spontaneous ventilation of patients during short anaesthesia procedures.

## CONTRAINDICATIONS

- Complete obstruction of the upper airway.
- In elective use: patients who have not fasted as well as patients with retained gastric contents.

No other contraindications are known.

## SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel with adequate skills in the use of the product.
- During MRI treatment, the pilot balloon of the tube can cause artefacts around the area where it is positioned. A non-clinical test verified that the tube can be used in an MRI environment.
- The product is suitable for conducting oxygen, air and anaesthesia gases. The anaesthesia gases used must not interact with the product materials (see material data).
- Laser or electrocaustic equipment may cause the product to be flammable.
- The tube and the syringe are for single use and must not be reprocessed. Reprocessing leads to impairment of the function of the products. Any reuse entails a potential risk of cross-contamination.
- The tube is sterile (ethylene oxide). The product must not be used if the packaging is damaged or if shelf life has elapsed.

## USE

- ▶ Choose the correct tube size (see table).
- ▶ Remove the valve opener.

### Perform the following visual and functional check:

- ▶ Check the tube for damage and loose particles.
- ▶ Check the inside of the tube for sufficient patency.
- ▶ Inflate the cuffs and check for leaks.

Damaged products must not be used.

Have a spare tube ready and prepared for immediate use.

- ▶ Deflate the cuffs completely (picture 1).
- ▶ Lubricate tip and posterior wall of the tube with water-soluble lubricant.
- ▶ Pre-oxygenation is advisable in all cases.
- ▶ Head position

The ideal head position for insertion is the "sniffing position". However, the tube's short length and its shape also allow it to be inserted with the head in a neutral position. For children, it is recommended (due to the anatomy) to have the head slightly extended with a chin lift manoeuvre.

### Insertion:

## CAUTION

Achieve an appropriate depth of anaesthesia before inserting the tube!

- ▶ Hold the tube like a pen in the area with the black lines (teeth marks) (picture 2).
- ▶ Open the mouth with the free hand (picture 3) and if possible, lift the chin. The tongue must not fold back during the downward movement of the tube.

**CAUTION**

Never use force to insert the tube! During insertion, the cuffs must not be damaged.

- ▶ Insert the tube either in the midline of the mouth or into the corner of the mouth while rotating the tube and advance to the upper teeth mark (picture 4).
- ▶ Use a syringe to inflate the cuffs with the recommended volume (see table) (picture 5).  
Thanks to the high-volume low-pressure cuffs, a seal can already be achieved with low cuff pressures (< 60 cmH<sub>2</sub>O).
- ▶ Connect the tube to a ventilation system (picture 6).

**CAUTION**

- ▶ Ventilation must be checked each time the patient is repositioned.
- Oxygen, air or anaesthesia gases can change the cuff volume or cuff pressure.  
Excessive cuff pressure can lead to swelling of the tongue.
- ▶ Reduce the cuff pressure to the Airway Leak Pressure and check regularly.

- ▶ Use an appropriate procedure to check ventilation and exclude gastric insufflation.
- ▶ If no sufficient ventilation can be verified, withdraw the tube with inflated cuffs between the teeth marks until ventilation is easy and free flowing (large tidal volume with minimal airway pressure).
- ▶ Check the cuff pressure again and adjust to the Airway Leak Pressure.
- ▶ Fix the tube to prevent dislocation. From size 2, the manufacturer offers a bite block with fixation tape.

The drain tube (see table) allows the insertion of a gastric tube or suction catheter.

**Removal:**

- ▶ Deflate the cuffs completely (picture 7).

Once the tube is in the correct position, it is well tolerated until the return of the protective reflexes.

- ▶ Remove the tube in an area where suction equipment and the ability for rapid intubation are present.

**PRODUCT SPECIFICATIONS**

| Tube size                          | #0                    | #1        | #2         | #2.5         | #3         | #4           | #5         |
|------------------------------------|-----------------------|-----------|------------|--------------|------------|--------------|------------|
| Weight/height                      | < 5 kg                | 5 - 12 kg | 12 - 25 kg | 125 - 150 cm | < 155 cm   | 155 - 180 cm | > 180 cm   |
| Colour coding ISO connector 15 mm  | transparent           | white     | green      | orange       | yellow     | red          | purple     |
| Recommended cuff volume            | 10 ml                 | 20 ml     | 35 ml      | 40 - 45 ml   | 50 - 60 ml | 70 - 80 ml   | 80 - 90 ml |
| Maximum cuff pressure              | 60 cmH <sub>2</sub> O |           |            |              |            |              |            |
| Drain tube                         | 10 Fr                 |           | 16 Fr      |              | 18 Fr      |              |            |
| Bronchoscopy via ventilation lumen | < 3.0 mm              |           | < 4.0 mm   |              | < 6.0 mm   |              |            |
| External diameter of the tube      | 9 mm                  |           | 14 mm      |              | 17.6 mm    |              |            |

| SCOPE OF DELIVERY                     |         |
|---------------------------------------|---------|
| Sets                                  | Syringe |
| Single set, child #0, #1              | 20 ml   |
| Single set, child #2, #2.5            | 60 ml   |
| Single set, adult #3                  | 60 ml   |
| Single set, adult #4, #5              | 100 ml  |
| Set of 10, child                      | ---     |
| Set of 10, adult                      | ---     |
| Emergency set, child #0, #1, #2, #2.5 | 60 ml   |
| Emergency set, adult #3, #4, #5       | 100 ml  |

| MATERIAL DATA                                    |  |
|--|--|
| Tube<br>Cuffs<br>Inflation line<br>Pilot Balloon | PVC (polyvinyl chloride)   |
| Tube connector                                   | PC (polycarbonate)   |
| Check valve and luer connection                  | PVC (polyvinyl chloride), PBT (polybutylene terephthalate), nitrile, stainless steel |
| Syringe  | PP (polypropylene), IR (isoprene)  |

### SHELF LIFE


Expiry date: see product label.

### STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

### DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international legal regulations.

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