



Instructions for use

Ambu[®] Redi-ACE Adult

Ambu[®] Redi-ACE Mini

MR

EN	This product is MR Safe
BG	Този продукт е безопасен при МР.
CS	Tento produkt je bezpečný při vyšetření magnetickou rezonancí.
DA	Denne produkt er MR-sikkert.
DE	Dieses Produkt ist MR-sicher.
EL	Αυτό το προϊόν είναι ασφαλές για MR.
ES	Este producto es compatible con RM.
ET	Toode on MR-kindel.
FI	Tämä tuote on MK-turvallinen.
FR	Ce produit ne présente aucun danger en cas d'IRM.
HR	Ovaj proizvod siguran je za upotrebu uz MR.
HU	Ez a termék MR-biztos.
IT	Questo prodotto è compatibile con RM.
JA	MR では安全です。
LT	Šis gaminys yra saugus naudoti su magnetinio rezonanso (MR) įranga.
LV	Šis ir MR drošs produkts.
NL	Dit product is MR-veilig.
NO	Denne produktet er MR-sikkert.
PL	Ten produkt można stosować bezpiecznie w badaniach rezonansu magnetycznego.
PT	Este produto é seguro para utilização em IRM.
RO	Acest produs este sigur pentru utilizare în IRM.
RU	Данное изделие безопасно для использования с МРТ.
SK	Tento produkt je bezpečný na použitie pri MR
SL	Ta izdelek je varen za uporabo z magnetno resonanco.
SV	Denna produkt är MR-saker.
TR	Bu ürün MR güvenlidir.
ZH	本产品可在 MR 中安全使用。

Content	Page
English	06
България	08
Česky	10
Dansk	12
Deutsch	14
Ελληνικά	17
Español	20
Eesti	23
Suomi	25
Français	27
Hrvatski	30
Magyar	32
Italiano	34
日本語	36
Lietuviškai	38
Latviski	40
Nederlands	42
Norsk	45
Polski	47
Português	50
Romania	52
Русский	55
Slovenčina	58
Slovenčina	60
Svenska	62
Türkçe	64
中文	66

Redi-ACE

Warnings

The Ambu Redi-ACE collar is intended for application by licensed health care practitioners as directed by a physician or other medical authority. Permanent injury may result from use other than as directed by a physician or other medical authority.

An extrication collar is only part of total patient immobilization.

It is important that practitioners are trained and practiced in the use and application of the device and that all instructions are read and understood prior to application of the collar.

All directions are to be considered as guidelines and not an attempt to define medical practice. Use the device only after review and approval from the supervising medical authority.

The collar is intended for Single Patient use only due to the risk of cross contamination.

Principles of Immobilization

Prevention and/or further extension of damage to the cervical spine and spinal cord is accomplished by maintaining the head and neck in a position of “neutral alignment.” The head and neck should not be flexed forward, extended backward, or cocked to one side or the other. The neutral alignment position prevents damaged vertebrae or foreign bodies from rubbing and/or cutting the spinal cord. Additionally, the neutral alignment position prevents the spinal cord from being twisted or compressed causing physical and/or ischemic injury to the spinal cord.

It is important that practitioners are trained and practiced in the use and application of the device and that all instructions are read and understood prior to application of the collar.

All directions are to be considered as guidelines and not an attempt to define medical practice. Use the device only after review and approval from the supervising medical authority.

Check the device before use.

Patient to be monitored while device is in use.

1. Intended use

The Ambu® Redi-ACE collar is a one-piece rigid cervical spine immobilization device designed to assist the rescuer with the maintenance of neutral alignment, prevention of lateral (side-to-side) sway and anterior-posterior (forward and backward) flexion and extension of the cervical spine during transport and routine patient care or movement.

2. Specifications

Storage temperature: Tested at -30 °C (-22 °F) and +70 °C (158 °F) according to EN 1789

Operational Temperature: -10 °C ~ +40 °C (14 °F ~ 104 °F)

3. Parts/Materials

Main panel:	Polypropylene	Velcro:	Polyamide
Sliding panel:	Polypropylene	Foam:	IXPE
Buttons:	Polyamide		

4. Sizing and Application

4.1 Sizing (1) (2)

- Measure the distance between an imaginary plane drawn horizontally and immediately below the patient’s chin and second horizontal plane drawn immediately on top of the patient’s shoulder.
- Compare this distance with the distance from the collar sizing line to the lower aspect of the plastic collar body (not the foam).

- Adjust the collar to the appropriate size.
- Simply pull the collar apart until the distance between the sizing line and the plastic collar body equals your finger measurement.
- Engage the safety locks by pushing DOWN the locking buttons.
- If the Ambu Redi-ACE needs to be resized, disengage the safety locks by lifting the locking buttons.
- The Ambu Redi-ACE may now be adjusted to the appropriate size.
- Now engage the safety locks by pushing DOWN the locking buttons.

4.2 Application (3)

- Proper application of any extrication collar requires two individuals. The first rescuer should maintain stabilization and neutral alignment of the head and neck.
- The second rescuer should slide the front of the collar along the patient’s chest and position the chin piece. The collar body should rest on the top of the patient’s shoulder and against the sternum without gaps. The patient’s chin must rest securely in the chin piece in neutral alignment position.
- While holding the front of the collar in place, wrap the back of the collar around the back of the patient’s head and neck and secure to the front of the collar with the Velcro® tab.
- Adjust the Velcro tab as necessary by holding the front cutout of the collar and tightening the Velcro tab.

4.3 Supine Application (4)

If the patient is supine, slide the back portion of the collar directly behind the patient’s neck.

The collar may be applied from either side. Once the back portion of the collar is evenly positioned behind the head, the front portion of the collar can be positioned under the chin. While maintaining the proper positioning of the collar with one hand, attach the Velcro with the other hand. Ensure that the loop Velcro tab is firmly attached and parallel to the fixed hook Velcro tab on the front of the collar.

4.4 Tube hook (5)

A nasal airway tubing can be attached to the hooks at both sides.

Warnings

Maintain neutral head and neck alignment during the entire procedure. Once the collar is properly applied and the patient is immobilized, the head should not be released. Further immobilize the head and neck with head blocks or a head immobilizer and secure the patient to the backboard with appropriate retention straps.

Caution

US federal law restricts this device to sale by or on the order of a physician.

Cleaning

To eliminate any risk of cross infection between patients it is recommended according to the local guidelines of medical waste to discard the product after use.

The surface of the product can be cleaned with water and a mild detergent or using a disinfecting agent. Rinse it thoroughly with clean water and allow it to dry before reuse.

Before reuse, make sure that the product is intact and fully operationable. Otherwise it should be discarded.

Warranty/Guarantee

Ambu warrants its products to be free from defects in workmanship and materials when purchased through Ambu or its authorized channels of distribution and used in accordance with recommended procedures. Should any product be found to be defective when unpacked, obtain a return authorization number from Ambu before sending it back. Repair or replacement will be made at no charge.

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Ambu is a registered trademark of Ambu A/S, Denmark.
Ambu A/S