Istruzioni d'uso
SHERPA/SHERPA PLUS – PORTABARELLA CAMPALI

IT

Instructions for use
SHERPA/SHERPA PLUS – FIELD STRETCHER CARRIERS

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EN

1. MODELS	8
2. INTENDED USE	8
2.1 INTENDED USE	8
2.2 TARGET PATIENTS	8
2.3 PATIENT SELECTION CRITERIA	8
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	8
2.5 USERS AND INSTALLERS	8
3. REFERENCE STANDARDS	8
4. INTRODUCTION	8
4.1 DEVICE LABELLING AND TRACEABILITY	8
4.2 SYMBOLS	8
4.3 WARRANTY AND SERVICE	8
5. WARNINGS/DANGERS	8
6. SPECIFIC WARNINGS	9
7. RESIDUAL RISK	9
8. TECHNICAL DATA AND COMPONENTS	9
9. PROPER USE	9
10. CLEANING AND MAINTENANCE	10
10.1 CLEANING	10
10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL	10
10.3 LIFE SPAN	10
11. TROUBLESHOOTING TABLE	11
12. ACCESSORIES	11
13. SPARE PARTS	11
14. DISPOSAL	11

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- SHERPA
- SHERPA PLUS

2. INTENDED USE

2.1 INTENDED USE

Sherpa and Sherpa Plus are support trolleys for field stretchers which can also be used for short trips within the equipped area. The Plus version is equipped with a system to keep the device in an inclined position. Compatible with series 100/200/300/Slim/Stock stretchers. Unsupervised parking is not permitted.

2.2 TARGET PATIENTS

The target patients are those for whom use of the field stretcher is foreseen.

2.3 PATIENT SELECTION CRITERIA

The selection criteria are those foreseen for use of the field stretcher.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the instructions for use.

2.5 USERS AND INSTALLERS

The intended users are those foreseen for use of field stretchers.
These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

The product must be used only by personnel trained in the use of this product and not on other similar products.

3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location.

4.2 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the instructions for use
	Manufacturer		Lot number
	Date of manufacture		Product code
	Unique Device Identifier		Caution: Federal law restricts the sale of this device by or on the order of a licensed professional (US market only)
			Production identification Alphanumeric code that identifies the production units of the device, composed of: (01)805771123 company prefix 000 progressive GS1 6 control number (11)200626 date of production (YYMMDD) (10)1234567890 lot number
	(01)8057711230006 (11) 200626 (10) 1234567890		

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase.**

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

Warranty and service conditions are available at <http://support.spencer.it>

5. WARNINGS/DANGERS

Product features

Use of the product for any purpose other than that described in the instructions for use is prohibited.

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to +50°C.
- Storage temperature: -10°C to +60°C.

General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

6. SPECIFIC WARNINGS


To use the product, you must also have read, understood and carefully follow all the instructions for use.

- ⚠ Never apply a weight greater than the maximum capacity of the device. When determining the applicable weight, the operator should consider the weight of the patient, device, equipment, accessories and any other object attached to the stretcher.
- ⚠ Perform simulations of use with a stretcher and a patient simulating load and accessories before putting the device into service..
- Never leave the patient on the device unsupervised, as they could get injured.
- Do not drag the device.
- Hold the stretcher firmly on the device if the patient is lying on it.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in these instructions for use, have been identified.

8. TECHNICAL DATA AND COMPONENTS


	DESCRIPTION
	1 Locking handle – Made of aluminium, this is used to attach the stretcher to the body
	2 Handle unlocking lever – Made of steel, used to unlock the locking handle
	3 Wheels – Made of polyurethane, they allow sliding on the ground
	4 Stands – Made of aluminium, these are folding ground support elements
	5 Tilting supports – Made of aluminium and only found on the Sherpa Plus; when extracted they allow one side of the device to be raised by a few centimetres.
	6 Central locking mechanism – Made of aluminium, allows the device to be locked in the open position
	7 Fastening strap – Made of nylon, it keeps the device closed when not in use

	SHERPA	SHERPA PLUS
Width when open	640 ± 10	640 ± 10
Width when closed	300 ± 10	300 ± 10
Length when open	995 ± 10	995 ± 10
Length when closed	670 ± 10	670 ± 10
Height	950 ± 10	950 ± 10
Materials	Al, Steel, Nylon	Al, Steel, Nylon
Weight	19 ± 0,5 kg	20 ± 0,5 kg
Load capacity	160 kg	160 kg


9. PROPER USE

Find a level area where you can place the device to ensure its stability.


To open the device:


Unhook the belt buckle that keeps the device closed, then spread the two halves of the frame apart.

Separate the two halves of the device completely, push the red sleeve in the direction slightly raising the joints, then release the sleeve when the joints are aligned on the same axis.

Open the support stand by turning it outwards.


Once the stand is fully open, keep the support stand (1) pressed to the ground with one foot, raise the device by about 5 cm by grasping the curved parts of the upper frame (2) then, with the other foot, press on the joint so as to lock the stand in the open position (3). Repeat the same manoeuvre on the opposite side of the device.


To close the device, grasp the curved part of the upper frame. Use one foot to push the joints of the stand described above in the opposite direction to the locking direction, then use the other hand to guide the stand into the closed position. Perform the same manoeuvre on the opposite side.
 Pull the red sleeve and push the joints down.
 Grasp the upper part of the frames and screw them together.
 Fasten the belt again.
 Failure to follow the closing sequence could lead to blocking of the device.
 During closing manoeuvres, make sure that the hands of other people who may be helping to keep the device stable are not within the handling area as they could be injured.

Stretcher positioning


Set the device in parking mode as described above.


Pull the locking handles out of their seats while simultaneously pressing the red button.

Place the stretcher at the centre of the device so that the unit is well balanced, taking into account the weight of the patient. Push the locking handles against the stretcher frame. This will lock the stretcher.

The Trendelenburg movement allows the stretcher to be tilted, thus raising the patient's lower limbs. Rotate the mechanisms on the stands approximately 180° so that they are on the other side of the rest position. Lower the stands to the floor and place the stretcher with the patient's limbs at the highest point.

To avoid accidental falls, **DO NOT ALLOW THE PATIENT TO GET ON OR OFF the stretcher when it is on the frame.**

The device can be moved by means of its wheels. During handling, the stands must be raised. Handling must always be carried out by at least two operators capable of keeping the unit balanced.

IT
EN

DE

FR

ES

PT

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues. The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.** Do not use cleaning agents containing chlorine / sodium hypochlorite as corrosion of components may occur.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If **disinfecting**, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL

Routine maintenance is the responsibility of the user, who must carry out checks before and after each use and in any case at least every 3 months:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Fulfilment of the requirements of the instructions for use in section Warnings and Specific Warnings.
- Fulfilment of the requirements of the instructions for use in the section on Proper use.

No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 10 years from the date of purchase. Beyond this time, the device must be disposed of.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

11. TROUBLESHOOTING TABLE

PROBLEMA	CAUSE	REMEDY
The stands do not lock when the joints are pushed.	The joints are damaged or the device has not been opened properly.	Check that the red sleeve has been correctly positioned. If the problem persists, immediately remove the device from service and contact the manufacturer.
Damage to the structure	Improper use or exceeded the maximum capacity	Remove the device from service and contact customer service.
The stretcher cannot be attached to the device	Locking handles not fully open Something is interfering with the locking system.	Check that you have opened the locking handles correctly. Make sure that nothing is interfering with the locking mechanisms. If damaged, immediately remove the device from service and contact the manufacturer.

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

12. ACCESSORIES

ST00003C	Rail bar for stretcher carriers	ST00004C	Side bag for stretcher carrier
----------	---------------------------------	----------	--------------------------------

13. SPARE PARTS

ST00013	Black wheel \varnothing 500 mm for Sherpa	ST00014	Sherpa stand bumper
ST00015	Stretcher locking belt for Sherpa		

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

IT

EN

DE

FR

ES

PT

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.
The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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