next**RESCUE**GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE


SpeedM®

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Sterile hemostatic dressing with HNT-mineral for control of traumatic bleedings

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until bleeding stops.

5. Check that the bleeding has stopped. Then apply a pressure bandage. Seek medical attention immediately and show the package.
6. Removing the product: after a maximum of 24 hours carefully remove the dressing from the wound. Rinse the wound with sterile solution.

Contraindications:

SpeedM® is not indicated for use on the eyeball due to the risk of dehydration or pressure damage. The product is not indicated for use in or in direct contact with the central circulatory system. SpeedM® must not be used in patients with a known allergic reaction to clay minerals.

Warnings:

⚠ Attention!

1. The product is non-absorbable and must be removed from the wound prior to wound closure.
2. SpeedM® Emergency Hemostatic Dressing is not intended as a substitute for careful surgical technique, the use of sterile procedures, and the proper use of ligatures or other conventional hemostatic procedures.
3. SpeedM® should be removed from the wound after a maximum of 24 hours.

Precautions:

Sterility is not guaranteed if packaging is damaged or opened. Dispose of the product if the packaging is damaged. For single use only. Protect from heat. Do not resterilize.

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Storage conditions:

Keep dry. Protect from heat, including storage in direct sunlight or in direct contact with heat sources. Store at a temperature of 25 °C (77 °F) or less.

Sterility and shelf life:

The product is sterilized by gamma radiation and is intended for single use only. Do not resterilize. Do not reuse. Do not use if the sterile bag is damaged or opened. Do not use the product after the expiry date indicated on the packaging.

Notes:

If serious incidents or unexpected events occur in connection with SpeedM®, please inform the manufacturer at nonconformity@speedcarenineral.de or the responsible authority.

Disposal:

Ensure that the product is disposed properly in accordance with national regulations.

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en INSTRUCTION FOR USE

Description:

SpeedM® Emergency Hemostatic Dressing is a sterile, 7.5 x 370 cm, radiopaque (detectable with X-rays), Z-fold hemostatic dressing and is packaged for asepti­c removal. This product contains halloysite (mineral).

Intended use:

SpeedM® is an emergency hemostatic dressing intended for temporary treatment of heavily bleeding traumatic wounds. The product is suitable for any patient to whom the intended purpose can be applied. The product is not suitable for children under 6 years of age.

Indications:

SpeedM® Emergency Hemostatic Dressing is intended for temporary treatment of heavily bleeding wounds and traumatic injuries such as cuts, lacerations and abrasions. For use by trained first responders.

Application:

1. Identify injury, locate active bleedings. Remove the folded gauze from the packaging. Care for the heaviest bleeding first.
2. Insert one end of the gauze down to the deepest point of the wound.
3. Then stuff and fill the wound cavity tightly from the bottom until the skin surface is reached. Depending on the size of the wound, several dressings can be inserted into the wound cavity. As the dressing is perforated, it does not matter whether the printed or unprinted side comes into contact with the wound first
4. Apply direct, firm pressure for at least 3 minutes or

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
- Identify injury, locate active bleedings. Remove the folded gauze from the packaging. Care for the heaviest bleeding first.
- Insert one end of the gauze down to the deepest point of the wound.
- Then stuff and fill the wound cavity tightly from the bottom until the skin surface is reached. Depending on the size of the wound, several dressings can be inserted into the wound cavity. As the dressing is perforated, it does not matter whether the printed or unprinted side comes into contact with the wound first.
- Apply direct, firm pressure for at least 3 minutes or until bleeding stops.
- Check that the bleeding has stopped. Then apply a pressure bandage. Seek medical attention immediately and show the package.
- Removing the product: after a maximum of 24 hours carefully remove the dressing from the wound. Rinse the wound with sterile solution.

- de Trainingsvideo (en)**
en Training video (en)
fr Vidéo de formation (en)
nl Trainingsvideo (en)
it Video di formazione (en)
es Video de formación (en)
pt Video de formação (en)
pl Film instruktażowy (en)
sk Treningove video (en)