



# SAM CHEST SEAL TRAINING

**SAM**<sup>®</sup>  
MEDICAL

**SAM**<sup>®</sup>  
MEDICAL

# WHY WE NEED CHEST SEALS

## WHY WE NEED CHEST SEALS



**Open chest wounds, specifically Tension Pneumothorax are the leading cause of potentially survivable battlefield deaths.**

# WHAT IS AN OPEN CHEST WOUND?

## WHY WE NEED CHEST SEALS



**Open chest wounds, specifically Tension Pneumothorax, are the leading cause of potentially survivable battlefield deaths. A valved chest seal should cover and protect the open chest wound while not preventing potentially lethal pressure from escaping the pleural space.**

1. Eastridge BJ, Mabry RL, Seguin P, Cantrell J, Tops T, Uribe P, Mallett O, Zubko T, Oetjen-Gerdes L, Rasmussen TE, Butler FK, Kotwal RS, Holcomb JB, Wade C, Champion H, Lawnick M, Moores L, Blackbourne LH. Death on the battlefield (2001-2011): implications for the future of combat casualty care. *J Trauma Acute Care Surg.* 2012 Dec;73(6 Suppl 5):S431-7.

2. Holcomb JB, McMullin NR, Pearse L, et al. Causes of Death in U.S. Special Operations Forces in the Global War on Terrorism: 2001-2004. *Annals of Surgery.* 2007;245(6):986-991.

A SWAT officer in tactical gear is shown applying a SAM Chest Seal to a person lying on the ground. The officer's uniform has a "SWAT" patch. The seal is a circular, adhesive device with a central white cap and a textured surface. The scene is dimly lit, suggesting an indoor or nighttime environment. In the background, a piece of equipment is labeled "E-100 PERFORMANCE FOC GENERATOR".

# SAM CHEST SEAL

<https://www.youtube.com/watch?v=H4RRmzjrsB0>



A man is lying in a hospital bed, appearing to be in a critical condition. He has a SAM medical device (a large, rectangular, multi-ported pad) applied to his chest. The device is secured with straps and has several tubes connected to it. A person's hands are visible, one resting on the man's shoulder and another near the device. The background shows the metal frame of the hospital bed and a blurred medical professional. The overall scene is dimly lit, emphasizing the medical emergency.

**ENGINEERED TO PRESERVE LIFE**

# ENGINEERED TO PRESERVE LIFE



# ENGINEERED TO PRESERVE LIFE



## TRUFLOW™ VALVE\*

Is composed of a reinforced dome which protects the chest wound and the internal one-way valve from external forces allowing air to flow from the chest cavity and preventing air ingress into the chest cavity.

\*Does not come with non-valved version

[https://www.youtube.com/watch?time\\_continue=2&v=IV17aBtQslg](https://www.youtube.com/watch?time_continue=2&v=IV17aBtQslg)

# ENGINEERED TO PRESERVE LIFE

## FORTIFIED GEL

Strong Hydrogel tightly adheres to the chest wall in the presence of blood, hair, sweat or sand. And if necessary allows lifting and re-application while maintaining superior adherence.



REVERSE SIDE

# ENGINEERED TO PRESERVE LIFE

## THIN EDGE LAYER

The thin edge layer resists lifting or peeling, maintaining a formidable seal



# ENGINEERED TO PRESERVE LIFE

## NIGHT VISION OPTIMIZED

Optimized for night vision, the packaging, pouch, dressing backing, and tabs provide ideal contrast for visibility in normal light, low light, and night conditions (using a night vision device).



# ENGINEERED TO PRESERVE LIFE

## OVAL SEAL

Larger than other dressings, its oval shape maximizes the surface area of the seal, and its transparency, elasticity and conformability ensure effective adhesion.





# ENGINEERED TO PRESERVE LIFE



## DUAL APPLICATION TABS

Tabs facilitate the proper placement and removal of adhesive seal, with the ability to vent and reseal.

# ENGINEERED TO PRESERVE LIFE

## PEEL / TEAR PACKAGING

Quickly and easily open SAM<sup>®</sup> Chest Seal packaging by peeling from its tabs or tearing along the chevron-marked areas.



# ENGINEERED TO PRESERVE LIFE



## ELASTICITY

Flexible polyurethane film allows the SAM Chest Seal to stretch and contract with normal respiratory movements, greatly improving adherence to the chest wall.

# FOUR CONFIGURATIONS

# FOUR CONFIGURATIONS



**COMBO**



**VALVE 2.0**



**NON-VALVED**



**WITH VALVE**

# COMBO



# NON-VALVED



# VALVED 2.0





# WITH VALVE




A soldier in camouflage gear is lying on a stretcher. A SAM chest seal is applied to their chest. The scene is dimly lit, suggesting a night-time or low-visibility environment. The text "WHY SAM CHEST SEAL?" is overlaid in large, white, bold letters across the center of the image.

# WHY SAM CHEST SEAL?

## SAM CHEST SEAL

**SAM<sup>®</sup> Chest Seal is engineered to treat, seal, and reseal open chest wounds under most circumstances. After simple and quick application, SAM<sup>®</sup> Chest Seal sticks relentlessly – through extreme heat or cold, no matter the elements. Meets CoTCCC-preferred features.**

# WHY SAM CHEST SEAL?



**OCCCLUSION -RESISTANT  
TRUFLOW VALVE**

WHY SAM CHEST SEAL?

STERILE R

**STERILE SEAL**


Minimize further  
contamination of the wound

## WHY SAM CHEST SEAL?

### **FORTIFIED GEL**

Strong Hydrogel adheres in the presence of blood, hair, sweat and sand, or when submersed in water.

# WHY SAM CHEST SEAL?



## **LARGE OVAL SIZE**

Larger than other dressings, its oval shape maximizes the surface area of the seal, and its transparency, elasticity and conformability ensure effective adhesion.

## WHY SAM CHEST SEAL?



### **NIGHT VISION OPTIMIZED**

Optimized for night vision, the packaging, pouch, dressing backing, and tabs provide ideal contrast for visibility in normal light, low light, and night conditions (using a night vision device).



# WHY SAM CHEST SEAL?



## 6-YEAR SHELF LIFE

Industry-leading shelf life ensures SAM<sup>®</sup> Chest Seal is always ready when you need it.

# WHY SAM CHEST SEAL?



A person in military uniform is lying on a patterned rug. A SAM chest seal is being applied to their chest. The seal is a rectangular piece of material with a central opening and a strap. The person's hands are visible, and they are holding the seal. The background is dark and out of focus.

# SAM CHEST SEAL APPLICATION

# COMBO




[https://www.youtube.com/watch?v=y1log\\_cklmM](https://www.youtube.com/watch?v=y1log_cklmM)



## STEP 1


**Quickly and easily open the sterilized SAM<sup>®</sup> Chest Seal packaging by peeling from it tabs or tearing along the chevron-marked areas.**



## STEP 2

---

**Clean the wounded area.**




## STEP 3

---


**Next, grip the application tab and remove the clear liner.**





## STEP 4


**Center the TRUFLOW Valve over the wound and place the dressing adhesive side down.**



## STEP 5

---

**Press firmly across its entire surface to ensure adhesion.**



## STEP 1

---


**Locate the wound and  
remove the casualties  
clothes.**



## STEP 2

---

**Clean the wounded area.**



## STEP 3


---

**Next, grip the application tab and remove the clear liner.**



## STEP 4

**Center the dressing over the wound and place the adhesive side down before pressing firmly across its entire surface to ensure adhesion.**



A person is lying on a stretcher, appearing to be in a medical emergency. A medical device is attached to their chest, and a hand is resting on their abdomen. The scene is dimly lit, suggesting an indoor medical setting.

# SUMMARY



**SUMMARY**

**STERILE**

**SUPERIOR ADHESION**

**OCCLUSION-RESISTANT TRUFLOW VALVE**

**LARGE SIZE**

**NIGHT VISION OPTIMIZED**

**6 YEAR SHELF LIFE**



# FAQS

## FAQS



## Why is the SAM Chest Seal a large oval shape?

The SAM Chest Seal is oval shaped and larger than other dressings. This maximizes the surface area of the Seal to assure effectiveness of adhesion. If the seal displaces from chest, it cannot protect the wound or the the air flow from chest. The large oval shape also facilitates coverage of the chest wound and surrounding tissue.

## FAQS



# Why does the SAM Chest Seal come with and without a valve?

The SAM Chest Seal comes in both versions because not all chest injuries are treated the same. The “valved” SAM Chest Seal provides a one-way valve to allow air flow out of the chest thereby preventing the life threatening condition; tension pneumothorax.

The non-valved version provides a large occlusive dressing. Primarily used for chest wounds when medical professionals have the means to carefully monitor for and treat tension pneumothorax using more definitive means.



## FAQS

### **Can the SAM Chest Seal be used in the prehospital setting and austere conditions?**

Yes, the SAM Chest Seal has been designed to work in many different environments including extreme temperatures (hot or cold), wet (blood, sweat, rain, etc.), in the presence of body hair, and dirt ( sand, dust). The SAM Chest Seal has been tested to assure you, the rescuer that even in truly austere conditions the SAM Chest Seal will work effectively and efficiently each and every time.

# FAQS



## Why are there so many notches in the SAM Chest Seal packaging?

The notches at the four corners allow quick access to the SAM Chest Seal in less than optimal conditions. No more searching for the tab or being forced to cut the package open. The top flap opening allows for the SAM Chest Seal to be opened such that the packaging itself provides a sterile field. In the field or in the hospital, the packaging of the SAM Chest Seal allows the most effective access when seconds count.

# REGULATORY

**EC REP**


**CE** 0050



**Rx ONLY**

**STERILE R**

Emergo Europe, Prinsessegracht 20,  
2514 AP, The Hague, The Netherlands



SAM Medical Products®  
27350 SW 95th Ave, Ste 3038  
Wilsonville, OR 97070 USA

**MORE INFO**

**sammedical.com**



**THANK YOU**

**SAM**<sup>®</sup>  
MEDICAL

**MORE THAN SURVIVAL™**