diagnostics
with passion

Serafol® AB0

 sifin diagnostics gmbh

de l en

sıfın

Serafol[®] AB0

de Information zum Gebrauch durch Fachpersonal	4
en Information for use by trained personnel	10

 sifin diagnostics gmbh
Berliner Allee 317-321
13088 Berlin, Germany

Phone: +49 30 700 144-0
Telefax: +49 30 700 144-30
E-Mail: info@sifin.de
www.sifin.de

2 °C

de**Serafol® AB0**

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25 °C

en**Serafol® AB0**

Bedside-card for identity proof
ready to use twin-card

INFORMATION FOR USE BY TRAINED PERSONNEL

Test purpose

Serafol® AB0 is used for an identity proof immediately prior to a blood transfusion (bedside test). This manual bedside test is the qualitative confirmation of earlier ABO- and D-blood typing of the recipient and ensures the compatibility of blood types of the recipient and the blood to be transfused. Thus possible mismatches can be detected. The test may additionally be performed with capillary-, vein-, or stored donation blood (red cell concentrate segment). Please refer to: valid directives for the collection of blood and blood components and those for application of blood products (haemotherapy).

Test principle

Test card coated with dried specific monoclonal reagents.

The test principle is a hemagglutination test for the detection of the respective red cell antigens. ABO blood types are defined by the presence or absence of the antigens A and B on red cells. If such antigens are present on the red blood cells, they will be agglutinated by the corresponding antibodies (positive reaction).

Twin card/ Single card

Serafol cards can be used either for ABO identity tests of recipients and blood donation (twin card) or for tests of 2 recipients (single card).

Composition

Serafol® AB0: The single card contains 3 reaction fields: one coated with Anti-A (clone A003) and one coated with Anti-B (clone B005) as well as one field for the blood to be tested or auto control.

Preservative: <0.1 % NaN₃

The antibodies are derived from culture supernatants of stable cell lines and demonstrate the consistent specificity and reproducibility characteristic for monoclonal antibodies.

Shelf life

Under proper storage conditions Serafol cards vacuum wrapped in foil can be used until the expiry date printed both on the foil and the outer label. Serafol cards should not be taken from their foil-wrapper until immediately before use. The Serafol card may only be opened and removed from the foil wrapper immediately before use. After opening, the card must no longer be stored.

Additionally needed reagents and materials

Isotonic saline solution or water (EU drinking water quality)

Pipettes (drop volume: 40-50 µl)

Test material

Capillary-, vein-, or stored donation blood (red cell concentrate segment)

If samples are collected from indwelling catheters, infusion systems and the like, ensure that additives from these systems are not taken up with the collected sample.

Test procedure

Method 1


1. Open foil wrapper and take out twin card.
2. When using as single card, fold twin card along perforation and tear apart prior to use.
3. Fill in recipient and blood donation data.
4. Add drops of material to be tested on reaction fields.

Recipient

Add one drop of recipient blood (40-50 µl) to each reaction field of one row (see figs. 1 and 2). Add one drop isotonic saline solution or water to each reaction field.

Donor (segment)


Add one drop (40-50 µl) donor blood to each reaction field of one row. Then add one drop of isotonic saline solution or water to each reaction field.

Auto control

The field "Blood" in each row is used to perform an auto control. Add one drop of isotonic saline solution or water plus one drop recipient or donor blood to the field (see also instructions "recipient" or "donor").

5. Continue with test procedure immediately. Do not let blood dry.
6. Mix each field with a new mixing stick for approx. 30 seconds until the reagent is completely dissolved. Only spread the mixture within the reaction field (see fig. 3). Use a new mixing stick for each reaction field in order to avoid antibody spreading. Use this procedure also for the auto control.
7. Keep the card in circular motion for approx. 30 to 60 seconds so that the drops circle in the reaction field (see fig. 4).

8. Observe for agglutination on the test fields and write test result on card.

Test results:

- Agglutination = positive reaction
- No agglutination = negative reaction

The four possible blood type combinations (A, B, AB and O) are listed in the following chart:

Anti-A	Anti-B	Auto control	Blood group
+	-	-	A
-	+	-	B
+	+	-	AB
-	-	-	O
+	+	+	Test void

If the Serafol card result does not match the recipient or blood donation **DO NOT TRANSFUSE**, immediately find the cause and contact a transfusion medicine specialist! To rely on the test result, a proper auto control is required. If the auto control field shows agglutination, the test result is not valid and must be re-checked.

Method 2

1. to 3. see **method 1**
4. Add 1 drop of isotonic saline solution or water to the reaction fields.
5. Add 1 large drop (70-80 µl) recipient or donor blood to one of the fields labelled "blood" on the card.
6. Dip a mixing stick as flat as possible into the blood drop and turn the stick until a ball of blood adheres to the end of the stick .
7. Dip the stick end with the blood into the drop of isotonic saline solution or water on the reaction field with Anti-A and mix until there agent is completely dissolved (approx. 15 seconds). Only spread the mixture within the reaction field. Repeat this procedure with a new mixing stick on the reaction field Anti-B in the same row.
8. Set card aside for approx. 30 seconds
9. Cautiously wave the card or turn to all sides. After 60 to 90 seconds observe for agglutination and write results on card.

Documentation

- The test result is to be recorded by the administering physician.
- For filing and documentation, the reaction mixes can be dried and covered by a self-adhesive, see-through foil cover. During the drying the cards should be stored horizontally and should not be moved.

Precautions and Warnings

- Manual techniques are to be performed according to the manufacturer's instructions. Each deviation from these instructions is the sole responsibility of the user.
- Used test materials must be discarded as hazardous material.
- The bedside test is no proper substitution for complete blood typing or serological compatibility tests (cross matching).
- Because of the reagents' biotechnological production the risk of contamination with infectious agents is almost impossible. As some content derives from animal material (fetal calf serum, stabilizer) all test reagents must be regarded as potential transmitters of infectious diseases and be handled accordingly.
- As the test reagents contain sodium azide avoid any contact with the skin or the mucosa!

Performance characteristics and limits of the procedure

The performance evaluation according to the Commission Implementing Regulation (EU) 2022/1107 laying down common specifications (CS) yielded the following calculated diagnostic values:

	positive	false negative	Sensitivity	negative	false-positive	Specificity
Anti-A	875 (343)	0	100%	960 (315)	0	100%
Anti-B	453 (139)	0	100%	1382 (519)	0	100%

The values in brackets involve provisions administered at the patient's bed.

Sources of errors and limitations

For cancer patients and those with haematopoietic stem cells disorders, changed reactivities may arise (see pathophysiology and therapy). Deviations must be clarified in the immunohaematology laboratory as a matter of principle.

Causes for false positive results

- Blood samples may occasionally react by rouleaux formation¹, which can be mistaken for a weak agglutination and may incorrectly be read as a positive result. This phenomenon has non-immunological causes and usually occurs when blood is directly and without isotonic saline solution or water added to the reaction fields and mixed. Rouleaux formation also occurs in heparin blood, blood from patients treated with plasma expanders (e.g. dextran or hydroxyethyl starch) as well as in blood from patients with plasmacytoma (high protein content, changes in protein content), oncological disease (pathological blood count) or coagulation dysfunctions.

For testing these patients, first 1 drop of isotonic saline solution and then a drop of blood of the same size should be added to the reaction fields. If rouleaux formation is suspected as the cause for false positive results or other reaction results for either the recipient or donor (auto control is not clear), the test should be repeated with method 2, which usually avoids this phenomenon.

- Drying out may be incorrectly interpreted as a positive reaction.
- **Perform test without delays.**
- Cold-reacting autoantibodies with high temperature amplitude may also cause a false positive result. These antibodies are detected during laboratory testing and have to be taken into account.

¹ Reference: Issitt, Peter, D.: Applied Blood Group Serology, 4th Edition 1998, page 1134

Causes for false negative results

- Blood drop too small or too large.
- **Repeat test precisely following instructions.**
- Expiry date has passed or card was stored improperly (over 30 °C for a longer time period). **Repeat test with new card.**
- Test was performed with red blood cell suspension.
- **Repeat test with vein-, capillary-, or stored donation blood** (red cell concentrate segment).
- Erythrocytes with weakly developed blood grouping (e.g. A_y, A_x) or with low antigen density are not always detected.
- Blood with a hematocrit under 15 % may cause weak or false negative reactions.

Vigilance

All serious incidents occurring in connection with the product must be reported to sifin diagnostics gmbh as manufacturer and the responsible authority of the member state in which the user and/or the patient is based.

Literature

- (1) Petershofen EK 2011. Über AB0-Blutgruppen, Transfusions-schemata und die Bedeutung von Bedside-Tests. Hämotherapie 17: 18-25
- (2) Giebel, Picker, Gathof. Evaluation of Four Bedside Test systems for Card Performance, Handling and Safety. Transfus med Hemother 2008; 35: 33-36

Package size

Serafol® AB0	Single use kit: 1 twin card, 2 cover foils, 2 mixing sticks	Catalogue no.: 803100, BG1721
--------------	--	----------------------------------

Symbol Explanation

LOT	Batch Code (Lot)		Use by YYYY-MM (MM = end of month)
REF	Catalogue number		Temperature limitation
IVD	In Vitro Diagnostic Medical Device		Consult instructions for use
 ABO	ABO system	BEDSIDE	Use by specialist medical personnel at the patient's bed
CLONE	Clone designation		Do not re-use
			Do not use if damaged package

Revision 1, Date of revision: 17/07/2023

 sifin diagnostics gmbh

Berliner Allee 317-321

13088 Berlin, Germany

Phone: +49 30 700 144-0

Telefax: +49 30 700 144-30

E-Mail: info@sifin.de

www.sifin.de

sifin