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User Manual
SPENCER KINETIX EN
Powered stretcher with platform

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EN

1.	MODELS	21
2.	INTENDED USE	21
3.	REFERENCE STANDARDS	21
4.	INTRODUCTION	21
5.	WARNINGS	22
6.	SPECIFIC WARNINGS	23
7.	RESIDUAL RISK	25
8.	TECHNICAL DATA AND COMPONENTS	25
9.	COMMISSIONING	26
10.	OPERATING CHARACTERISTICS	27
	10.1 Drilling	27
11.	PROPER USE	29
	11.1 Stretcher Operating commands	29
	11.2 Platform operating commands	30
	11.3 Loading and unloading the stretcher	32
	11.4 Control panel	33
	11.5 Errors key of control panel	33
	11.6 Breakdowns	34
12.	CLEANING AND MAINTENANCE	34
	12.1 CLEANING	35
	12.2 ROUTINE MAINTENANCE	35
	12.3 PERIODIC OVERHAUL	36
	12.4 SPECIAL MAINTENANCE	36
	12.5 LIFE SPAN	36
13.	TROUBLESHOOTING TABLE	36
14.	SPARE PARTS	37
15.	DISPOSAL	37

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- SPENCER KINETIX – POWERED STRETCHER WITH PLATFORM

Devices above need to be used mandatorily together.

2. INTENDED USE

Ambulance stretchers are the primary means of transporting sick and/or injured people in transit safely and comfortably in a lying position. It is not foreseen that the patient be able to intervene on the device, nor is it intended for prolonged storage or use as a hospital bed. The stretcher must be used with the dedicated fixing systems. It is not foreseen that the patient be able to intervene on the device.

TARGET PATIENTS

There are no particular indications related to the patient group.

The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity of the device. If paediatric subjects must be transported, it will be in the role of the rescuer to determine whether the strap systems are suitable for immobilization or if it will be necessary to use another device.

PATIENT SELECTION CRITERIA

The intended patients are those for whom ambulance transport is necessary.

CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

USERS AND INSTALLERS

The intended users are persons trained in first aid procedures and in use of medical equipment in an emergency medical service (EMS) environment.

Possible users also include emergency vehicle outfitters who may use the product before commissioning or during maintenance of the vehicle on which the stretcher is used.

User training

Note : despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment may sometimes differ significantly.

The best instructions are the continuous practice of use under the supervision of competent and trained personnel.

- Regardless of your level of experience with similar devices in the past, you should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.
- The product must be used only by personnel trained in the use of this product and not on other similar products.
- The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

Note : Spencer Italia S.r.l. is always available for training courses. Document user training using the form in Annex A of this manual or other appropriate forms.

Installer training

Device installation must be carried out by qualified personnel who have been trained and certified by Spencer Italia S.r.l.

The installer must strictly adhere to these instructions as well as to the state of the art for vehicle installations.

3. REFERENCE STANDARDS

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

REFERENCE	DOCUMENT TITLE
REGULATION (EU) 2017/745	REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices.
UNI EN ISO 1865-1	Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
UNI EN 1865-3	Patient handling equipment used in road ambulances - Part 3: Heavy duty stretcher
UNI EN 1789	Medical vehicles and their equipment - Road ambulances

4. INTRODUCTION

4.1 USING THE MANUAL

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.

Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.

Spencer products User Manuals can be downloaded from the site www.spencer.it/en/emergency-medical-solutions/user-manuals or by contacting the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

IT

EN

DE

FR

4.2 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

If the assigned Lot/SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from Spencer Italia S.r.l., please register the device at service.spencer.it, or inform Customer Service (see § 4.4).

4.3 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Medical device		See the user manual
	Manufacturer		Product code
	Date of manufacture		Serial Number
	Unique Device Identifier		

	Production identification Alphanumeric code that identifies the production units of the device, composed of: (01)805771123 company prefix 000 progressive GS1 6 control number (11)200626 date of production (YYMMDD) (21)1234567890 SN		
(01)8057711230006 (11) 200626 (21) 1234567890			
	Waste symbol for electrical and electronic equipment “WEEE” (Waste Electrical and Electronic Equipment). Recycle: Electronic equipment. DO NOT THROW GARBAGE		“Dangerous voltage” symbol. To warn of electricity.
	“Double insulation” symbol. The protection of appliances marked with this symbol is ensured by double insulation and does not require a safety electrical connection to earth (ground).		Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)
	Warning: crushing of hands		Indicates the temperature limits; The upper and lower temperature limits should be indicated adjacent to the lower horizontal line

4.4 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.
For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541111, fax +39 0521 541222, e-mail service@spencer.it.

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

Warranty and service conditions are available at www.spencer.it/en/support/terms-and-conditions.

Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal www.spencer.it/en/services/after-sales, available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

5. WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website www.spencer.it on the specific product page.

Product features

- Use of the product for any purpose other than that described in the User Manual is prohibited.**
- Before each use, always check the conditions of the product, as specified in the User Manual. In the event of faults/damage that could compromise its functionality/safety, immediately remove it from service and contact the Manufacturer.
 - If the product is found to be malfunctioning, immediately use a similar device to ensure continuity of ongoing operations.
 - The product must not be tampered with (modification, tweaking, additions, repair). Should these operations be performed, we decline any responsibility for incorrect operation or any damage caused by the product itself; moreover, the CE certification (when required by law) and the product warranty shall be null and void.
 - Anyone who modifies or causes products manufactured by Spencer Italia S.r.l. to be modified or re-processed in such a way that they no longer serve the intended purpose or no longer provide the intended performance shall comply with the conditions valid for the first placing on the market.
 - When using the devices, position and adjust them in such a way that they do not hinder operator works or the use of any other equipment.
 - Be sure to take every precaution to avoid hazards from contact with blood or body secretions, if applicable.
 - Always comply with the maximum capacity indicated in the User Manual. Maximum load capacity means the total weight distributed according to human anatomy. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product.
 - Make sure that operators are in good physical condition before lifting, as listed in the User Manual.
 - The maximum weight, which weighs on each operator, must comply with local health and safety requirements.**
 - Avoid contact with sharp objects.
 - Device installation must be carried out by qualified personnel trained and enabled by Spencer Italia S.r.l. The times and methods of such shall be agreed upon between the customer and our Sales Offices.
 - Operating temperature: from -10°C to + 50°C.

Storage

- The product must not be exposed or come into contact with thermal sources of combustion or flammable agents, but must instead be stored in a dry, cool place, away from light and sun.
- Do not store the product under other more or less heavy materials that may damage the structure of the product.
- Store and transport the product with its original packaging, otherwise the warranty shall be invalidated.
- Storage temperature: -20°C to +60°C.

Regulatory requirements

As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

- Promptly and in detail notify Spencer Italia S.r.l. (already in the quotation request phase) about possible fulfillments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act with due care and diligence to help ensure compliance with the general safety requirements of the devices placed on the market, providing end-users with all the information necessary to carry out periodic revisions on the supplied devices, exactly as indicated in the User Manual.
- Participate in safety checks on products**, placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.
- Without prejudice to the above, the Distributor or End-User shall assume wider liability related to non-compliance with non-fulfilment of the above-mentioned obligations, with consequent obligation to indemnify and/or hold Spencer Italia S.r.l. harmless from any possible injurious effect.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

General warnings for medical devices

- The user must carefully read the following in addition to the general warnings.
- The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point.
 - Qualified personnel and at least two operators must be present during use of the device.
 - Follow the internal procedures and protocols approved by your organisation.
 - Disinfection operations must be carried out in accordance with the validated cycle parameters, as stated in the specific technical standards.

6. SPECIFIC WARNINGS

- To use the product, you must also have read, understood and carefully follow all the instructions in the user manual.
- Follow approved Emergency Medical Service procedures for patient immobilization and transportation.
 - Follow approved Emergency Medical Service procedures for patient positioning and transportation.
 - Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
 - Before any type of handling, make sure that operators have a secure grip on the device.
 - Do not drag the device on uneven surfaces.
 - Do not lift with a crane or other mechanical lift.
 - Do not use drying machines.
 - The device is a transport device and cannot be used as a stationing device.
 - Do not use with devices other than those expressly approved by the manufacturer.
 - Practice with a device without a patient to make sure you are familiar with the manoeuvres.
 - For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).
 - Make sure that the patient is adequately immobilized before loading him/her on the stretcher. A failure to immobilize can cause serious damage to the patient.
 - Make sure that the sheet does not interfere with any stretcher handling or control mechanisms.
 - Do not move the device if the weight is not properly distributed.
 - Always use straps that are anchored to the stretcher frame to ensure patient safety.
 - Use only the perimeter frame or the handles to move the stretcher and not the sidebars, beds or other points not intended for this purpose.
 - Avoid excessive force when loading the stretcher into the ambulance: unnecessary force can cause damage and can adversely affect the fastening system.
 - Avoid excessive force when loading the stretcher into the medical vehicle: unnecessary force can cause damage and can adversely affect the functionality of the stretcher itself.
 - Hold the device firmly if the patient is lying on it.**
 - Parking brakes are aids for the operator but they do not in any way replace his or her supervision.**
 - Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. If you cannot clear the way, choose an alternative route.
 - Condensation, water, ice and accumulations of dust can affect correct operation of the device, making it unpredictable and causing a sudden change in the weight the operators have to support.
 - For height differences greater than 10 mm, the device must be lifted, taking care to grasp it from the structure and not from the sidebars/beds or other points not intended for this purpose.
 - Once the wheels of the loading trolley have been placed on the ambulance floor, the wheels of the front leg must have a ground clearance between 5 and 10 cm, which allows the front leg to be opened and locked safely. Check the height of the ambulance bed after each use; if it is altered, the trolley must be adjusted immediately by the manufacturer or by a specialised technician authorised by the manufacturer. If not, the manufacturer cannot be held responsible for correct functioning of the device or any damage caused by it.
 - If the vehicle is equipped with pneumatic or hydraulic suspensions, the loading height must be adjusted taking into account the worst-case conditions and/or the working conditions foreseen by the outfitter.
 - Any problems regarding the use and/or safety risks associated with this system are not the responsibility of the manufacturer.
 - Inadequate installation of the loading floor can lead to buckling and subsequent damage to the welds of the front legs.**
 - Improper installation of the loading floor can lead to abnormal operation of the device and cause harm to the patient and user.**
 - Do not alter or modify the stretcher to adapt it to the emergency vehicle: doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the Manufacturer from all liability.
 - The product can only comply with regulation EN 1789 if used with the dedicated fixing system. The use of fasteners not approved by the Manufacturer is therefore prohibited. Non-approved fixing systems may alter the structural and functional characteristics of the device.
 - Do not operate the variable height system without carefully assessing the weight of the stretcher with the patient and any accessories. Operators should be able to fully support the load when changing heights. Incorrect judgement may result in the stretcher falling unexpectedly and posing a risk to the patient and operators.
 - Do not place any magnetic parts between the stretcher and the fixing system as they may interfere with the stretcher's attachment and release systems.

- Do not place limbs and/or objects between the legs and the frame, in the vicinity of the leg lifting pistons or in general between moving parts as this may result in crushing injuries.
- If the twist system has been activated, make sure that the wheels are locked again before loading the stretcher into the ambulance as this can be risky and difficult with the front wheels unlocked.
- Movement with four swivel wheels can be very difficult on sloping or uneven ground. Carefully assess the conditions of use before unlocking the front wheel pivot.
- If the device is used in conjunction with restraint systems such as spine boards and/or vacuum mattresses, make sure that the patient is secured to the stretcher and the restraint system so that the patient is safe while the vehicle is moving. In case of doubt about the procedures to be used, refer to the operating protocols of your local first response center.
- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Always set the stretcher load height to the proper stop height before operation.
- Always have a certified mechanic, familiar with ambulance vehicle construction, install the stretcher support. Consult the vehicle manufacturer before install the stretcher support. Make sure that the installation does not damage or interfere with the brake lines, oxygen lines, fuel lines, fuel tank or electrical wiring of the vehicle.
- Always charge the batteries before using the product. An uncharged or depleted battery may cause poor product performance.
- Only charge the stretcher's batteries with the charger supplied by Spencer, otherwise it could damage the batteries and will lose their warranty.
- Always clear obstacles that may interfere and cause injury to the operator or patient before operating the stretcher.
- Always remove the battery if the stretcher is not going to be used for an extended period of time (more than 24 hours)
- It is important to maintain a minimum of charge cycles each month to maintain the good capacity of the batteries.
- Do not leave the batteries uncharged for more than 2 MONTHS. They could stop working.
- Do not steam clean or ultrasonically clean the product.
- Do not clean, service, or perform maintenance while the product is in use.
- Always use authorized parts to avoid the risk of product damage.
- Always check hoses and lines regularly to avoid damage to the stretcher. Check and tighten loose connections. Hydraulic lines, hoses, and connections can fail or loosen due to physical damage, kinks, age and environment exposure
- The use of accessories, transducers, and cables other than those specified or provided by the manufacturer, could result in increased electromagnetic emission or decreased electromagnetic immunity and result in improper operation
- In repeatedly raising and lowering maneuvers of up to 20 uninterrupted cycles, we must make a mandatory stop for at least 15 seconds so that the electro-hydraulic system of the stretcher recovers.
- To load the stretcher, the rail must always be fully extracted outside the ambulance.
- The electronic system of the stretcher will automatically switch off if no operation has been carried out within 15 minutes.
- When the stretcher is fully loaded on the rail, it will stop after 3 seconds.
- Use a protective coating for the stretcher mattress.
- The front wheels of the stretcher must remain locked when loading / unloading the stretcher. If not, important elements of the stretcher may be damaged.
- Always keep your hands clear of the moving parts when you load, unload or change the height position of the stretcher
- Always use both hands when transport the stretcher
- Do not attempt to operate the stretcher when is loaded on the stretcher support
- Do not remove the battery when the stretcher is active
- Always operate the stretcher when all persons are clear of the mechanisms. Entanglement in electric stretcher mechanisms can cause serious injury
- Do not allow untrained assistants to assist in the operation of the product
- Do not ride on the stretcher support
- Do not transport the stretcher sideways to avoid the risk of tipping. Always transport the stretcher in a lowered position, head end or foot end first to minimize the risk of tipping
- Always conduct the patient monitoring when the stretcher is idle- If you hydraulically raise or lower the product you may temporarily affect electronic patient monitoring equipment
- Do not leave the patient unattended. Hold the product while a patient is on the product
- Do not apply wheel lock when a patient is on the product or when you move the product to avoid the risk of tipping
- Do not use handrails as a patient restraint device
- Always transport the stretcher at a lower height to reduce the risk of a stretcher tip. If possible obtain additional assistance or take alternate route
- Always avoid high obstacles, such as curbing, steps, or rough terrain to avoid the risk of the product tip-ping over
- Always use all the belts to secure the patient on the stretcher. An unfastened patient may fall from the stretcher
- Do not attempt to open the battery pack for any reason, to avoid the risk of electric shock. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling
- Always avoid direct contact with a wet battery or battery enclosures. Contact may cause injury to the patient or .
- Always wipe the product with clean water and dry after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause a premature corrosion of critical components
- Always relieve pressure before you disconnect hydraulic or other lines. Escaping fluid under pressure can penetrate the skin and cause serious injury. Tighten all connections before you apply pressure. If an accident occurs, see a doctor immediately
- Do not use bare hand to check for hydraulic leaks
- The electrical installation of the rail must be carried out with the ambulance electrical systems disconnected.
- Do not lubricate the rail and / or stretcher when they are connected to the electrical network.

6.1 PHYSICAL REQUIREMENTS FOR OPERATORS

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands. **Operators must be able to provide the necessary patient care. Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device. For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).** The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.


7. RESIDUAL RISK


The residual risks listed below have been identified solely with reference to the intended use of the device.

- Use by untrained personnel may result in injury to the patient, rescuer or third parties.
- Inadequate disinfection procedures may lead to risk of cross-infection.
- Partial opening of the legs may result in the device falling to the ground. Make sure that the legs are properly locked before any handling and that the pistons are fully extended and stable.
- If the stretcher is not locked on the fixing system or if it is not correctly positioned, this can result in dangerous movements, especially in the event of strong medical vehicle deceleration, which can cause damage to the patient and the operators. Always make sure that the locking system is properly engaged.
- Failure to observe the warnings for operators can lead to risks of crushing caused by the moving mechanisms.
- Accidental operation of the variable height system may result in the stretcher falling, causing injury to the patient and/or operators. Make sure that the unlocking handle is not accidentally activated.
- The operators should prepare to fully load the stretcher with the patient and accessories before operating the variable height unlocking control. Operating this control without applying sufficient support force will result in the stretcher falling unexpectedly and causing injury to the patient or operators.
- **Failure to read and understand the product instructions may result in injury to the patient and operators.**




8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without notice.

Feature	Dimension
Length (mm)	2000
Width (mm)	590
Height (mm)	Min. 420 Max. 1300
Backrest angle	75°
Fowler angle	17°
Trendelenburg angle	30°
Max Loading height (mm)	1095
Input power supply	12 VDC - 72W
Power consumption	768 W
Load capacity (kg)	300
Weight without accessories fitted as standard (kg)	88
IP Classification	IPX6

Feature	Dimension
Length (mm)	2080
Width (mm)	600
Front movement (mm)	1640
Sideways movement (mm)	300
Weight without accessories fitted as standard (kg)	140 kg
Loading capacity (kg)	300 kg

9. COMMISSIONING

For first use, check that:

- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- The medical vehicle is equipped with the Spencer Kinetix platform dedicated to the stretcher
- The stretcher floor is properly levelled
- The stretcher floor is wide and long enough to accommodate the stretcher and its accessories without becoming an obstacle
- The fastenings must keep the stretcher attached to the structure of the vehicle

Do not modify the stretcher structural, leverage or traction parts in any way as this could cause damage to the patient and/or the rescuers.

Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

To facilitate insertion of the stretcher into the ambulance, eliminate sharp edges on the edge of the ambulance loading floor. The stretcher should be secured using Spencer hooks in such a way as to prevent any movement during transport in the ambulance, even in difficult driving conditions. Practice with a stretcher without a patient before putting the device into service. **For subsequent use, perform the operations specified in paragraph 12.**


If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and will void the warranty and release the Manufacturer from all liability.

10. OPERATING CHARACTERISTICS

10.1 DRILLING

Before putting the stretcher support into operation, ensure that each and every anchorage points are placed and clamped.

IMPORTANT FOR INSTALLATION

- Screws to be attached on the floor should be M12 according to norm DIN-7991 (8 units)
- Before drilling the ambulance floor, verify that the lateral movement of the stretcher support is placed at 0 position
- Install it upon the limits indicated on the drawing and make sure to fix the set on the rein-forced chassis are
- The electric installation of the set must be completed with the electrical system of the ambulance disconnected

10.2 POWER CONNECTION

To make the connection follow the steps below.

Be particularly careful when the connecting the stretcher support to the required voltage of 12V and do not leave cables in passageways and or moving elements. It can provoke breakages, cuts and system failures.

To make the connection locate the black cable placed on the front side of the support stretcher.

Remove the plastic end of the BROWN (+) y BLUE (-) cables with a thread peeler to make the connection.

Use a power strip and a flathead screwdriver to connect the BROWN (+) y BLUE (-) wires together with the outlet of the ambulance.

Make sure you have the power line properly protected and ensure the power along the line is stable.

10.3 START UP

To start using the electric stretcher it is necessary to charge the batteries. To do this, locate the external charger connector on the right side of the battery cabinet, on the back of the stretcher. Using the charger supplied by Spencer, charge the batteries to their maximum charge level (approximately 30 minutes).

Charge only the batteries of the stretcher with the charger supplied by Spencer, otherwise the batteries could be damaged and lose their guarantee.

11. PROPER USE

11.1 STRETCHER OPERATING COMMANDS

1. Foldable and reclining handrails
2. Retractable belts on the backrest
3. Touch control panel in the middle
4. Batteries cabinet with lock
5. Four lockable wheels
6. Telescopic frame for maneuvering
7. Bed with fowler and trendelenburg
8. Folding and adjustable IV pole
9. Powerful LED lighting on the sides
10. Emergency strap
11. Belts with front locking push button
12. Automatic integrated charging system
13. Front wheel lock pedal
14. Adjustable headrest

Adjust the backrest inclination.
Have 6 positions. Pressure the red handle and accompany the backrest to desired position.

Adjust maximum high of the stretcher
To adjust the maximum high of the stretcher loosen the screw and spin to desired position.


Waterproof batteries cabinet
To open the batteries cabinet, unscrew the 4 screws

Four lockable wheels
Spin the red handle to block or unblock the turn of the stretcher's wheels

Telescopic frame
Place the backrest in the vertical position, trigger by lateral knob and pull out or close the frame

Bed with fowler and trendelenburg position
The bed set has two red buttons to adjust the position

Folding and adjustable IV pole included
Turn the knob to raise the support rod and turn the upper ring to open it.

Powerful LED lighting on the sides
Press the button LIGHTS from control panel for activate the lights in the blinking or fixed mode

Emergency strap for lowering the stretcher
If you had to fold the stretcher in the manual way pull the emergency strap to fold it

Belts with front locking push button
To unfasten the belts, push the red button. You must always fasten the patient

Automatic integrated charging system
When you keep the stretcher on the support, the batteries of the stretcher will start to charge


Connector for external charger
The stretcher has a connector on the side of the battery box to charge it with an external charger

Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)


11.2 PLATFORM OPERATING COMMANDS

1. Pull-out lever for removable rail
2. Third point stretcher lock
3. Stretcher loading point
4. Stretcher central hooks
5. Lever to pick up the removable rail
6. Ball bearing to support removable rail
7. Lever to collect the removable rail from inside of the ambulance
8. Side movement lever
9. Safety hooks to grip the stretcher
10. Scoop stretcher housing module (if available)

LATERAL DISPLACEMENT – Hand lever
Operate the lever located on the rear of the right side of the support, making a slight effort with foot and accompanying that movement by sliding with hands. The total lateral displacement is 300 mm.

LONGITUDINAL DISPLACEMENT – Hand lever
Operate the lever located on the front of the right side of the support by applying a slight effort with the foot and accompanying this movement by sliding your hands. The total longitudinal displacement is 1640 mm.

11.3 LOADING AND UNLOADING THE STRETCHER

⚠ To load the stretcher, the support of stretcher always must be placed totally pull out of the ambulance.
Starting from the position in which the support has been left after the loading of the stretcher into the ambulance (1):



UNLOADING

- Pull the lever of the anchor point towards the left (A).
- Automatically the stretcher will be free. Push it out the ambulance (B) until the support are totally extract (2).
- Use the control panel to power the stretcher. WAIT 2 SECONDS to activate the system and later push the UP button (C).
- Once the stretcher will be in its highest position, which will be free (4).
- Move away the stretcher in straight line until being completely out of the support.
- Again, pick the support up into the ambulance by the lever placed on right side (F).

LOADING


- Pull the support out the ambulance by lever of the anchor point (A).
- Push the support by lateral handles until its maximum.
- Move the stretcher closer towards the support in straight line (5) make sure the wheels are placed like the image (D).
- Push the DOWN button on the control panel of the stretcher (E) and after push right lateral lever of the support (F) to put the set into the ambulance.
- Secure the stretcher is fasten on the anchor point (A) and anchor point is on close position.
- Verify the stretcher are charging, pushing the POWER button.

- **⚠** The front wheels of the stretcher must remain locked when loading / unloading the stretcher. If not, important elements of the stretcher may be damaged.
- When the stretcher is fully loaded on the rail, it will stop after 3 seconds.

11.4 CONTROL PANEL

All the electrical functions of the stretcher are controlled from the touch control panel placed on the back of the stretcher. See its functions in the following graphic.


BATTERY LED: the duration of the batteries is 4 months if the stretcher is stopped, 2 weeks if it is on but without maneuvering or 100 maneuvers with weight.
MAINTENANCE LED: the amber light will come on after 3,000 hours of use to perform maintenance on the stretcher. In that case, you should contact your distributor or Spencer.
If 4000 hours of operation are exceeded, the ERROR LED will light.

11.5 ERRORS KEY OF CONTROL PANEL

The led indicator of batteries status indicates possible errors of the system. In the following table, it is described the error according to its combination:

1. ERROR OF CONNECTION BETWEEN STRETCHER AND SUPPORT
2. ERROR FOR ELECTRIC SUPPLY ON THE SUPPORT
3. ERROR FOR OVERCHARGING BATTERIES
4. ERROR TIMEOUT FOR LOADING
5. ERROR MOVING UP FOR EXCESSIVE TIME
6. ERROR MOVING DOWN FOR EXCESSIVE TIME
7. ERROR IN CONTROL DRIVER OF STRETCHER
8. ERROR IN ENGINE TO EXCEED MAX. INTENSITY
9. ERROR IN ENGINE FOR LOW TENSION
10. ERROR FOR EXCEED LOAD

11.6 BREAKDOWNS

HYDRAULIC SYSTEM FAILURE

In case of failure of the hydraulic system, it has a strap that allows the stretcher to be folded manually. To do this, locate the strap on the top of the actuator.

Pull the strap and the stretcher will descend until its lower position.
Once the stretcher is folded, you should proceed to disassemble the safety hooks of the support as shown in the picture.
Once the stretcher is loaded into the support, you should place safety hooks in its initial position.

FAILING ONE OF THE BATTERIES

In case of failing one of the batteries, the stretcher will work perfectly.
The system is prepared to work with one battery. The second one will come in automatically when the control software detects that the first battery is about to run out.
However, if you have to change the batteries, you should proceed to open the compartment by unscrewing the four screws. Loosen the side screws that secure the battery holding cover and push it upward to remove the batteries.

 Batteries are delivered with a minimum charge. Before using the stretcher, you should previously charge the batteries with the external charger supplied with the stretcher to its maximum charge level. For a full charge, you should wait about 30 min.

12. CLEANING AND MAINTENANCE

Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised to do so; moreover, doing so will invalidate the warranty.

- The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.
- Establish a maintenance schedule, periodic inspections and extend the average life span, if foreseen by the Manufacturer in the User Manual, identifying a reference person who meets the basic requirements set forth in the User Manual.
- The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.
- Repairs of products manufactured by Spencer Italia S.r.l. must be carried out by the Manufacturer, who shall make use of specialised internal or external technicians who, using original spare parts, shall provide quality repair service in strict compliance with the technical specifications indicated by the Manufacturer. Spencer Italia S.r.l. declines any responsibility for any direct or indirect damage which is a consequence of improper use of spare parts and/or any repair work carried out by unauthorised parties.
- Use only original or Spencer Italia S.r.l. approved components/spare parts and/or accessories in order to carry out all operations without causing alterations or modifications to the product.
- All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested.
- Cleaning, provided for reusable products, must be carried out in accordance with the Manufacturer's instructions in the User Manual in order to avoid the risk of cross-infection due to the presence of body fluids and/or residues.
- The product and all its components must be washed and left to dry completely before storage.
- If the product requires lubrication, this must be done after cleaning and complete drying.
- The dedicated stretcher fixing system is equipped with a system that signals the need for stretcher maintenance and of the fixing itself linked to the cycles of use. Strictly observe the required maintenance intervals.
- Document maintenance using the form in Annex B of this manual or other appropriate forms.

12.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.
The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Any metal parts exposed to external agents undergo surface treatments and/or coating in order to obtain better resistance. Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers. Do not use cleaning agents containing sodium hypochlorite as corrosion of components may occur.**
Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.
If disinfecting, use products that do not have a solvent or corrosive action on materials constituting the device, in addition to being classified as medical-surgical devices. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

12.2 ROUTINE MAINTENANCE

Establish a maintenance program and periodic checks, identifying a designated reference person. The person entrusted with device maintenance must ensure the basic requirements envisaged by the manufacturer in the following paragraphs. All routine and special maintenance activities and all general overhauls must be recorded and documented with the relevant technical intervention reports. This documentation must be kept for at least 10 years after the end of the device's life and must be made available to the competent authorities and/or the Manufacturer when requested.

To ensure the traceability of the products and protect maintenance and service procedures on your devices, Spencer has made the SPENCER SERVICE portal www.spencer.it/en/services/after-sales, available to you. From this site, you can view the data of the products in your possession or placed on the market, monitor and update schedules for periodic checks and view and manage special maintenance.

Routine device maintenance must be carried out by operators with specific qualifications, training and instruction in the use and maintenance of the device.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking, maintenance and cleaning operations.


Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Correct tightening of all nuts, bolts and screws
- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- Welded areas are intact, without any cracks or breaks
- Moving parts, wheels, levers, handles are intact and work properly
- Lubrication of moving parts
- Conditions of use of wheels and brake system
- The wheels are correctly fixed, are stable and turn properly
- The wheels are free of debris
- The device opens and locks properly
- The device opens and closes properly
- Functioning of springs
- The stretcher easily enters the ambulance
- The medical vehicle is equipped with a Spencer fixing system intended for the stretcher
- The coupling between the fixing system and the stretcher is suitable to guarantee the safety of the fixing.

CONTROLS AND SELECTORS: Due to a heavy and continued use of the moving parts, such as controls and selectors, check them periodically for a proper functioning. There may be failures on micro-switches and cabling; examine the fixations and electric or mechanical connections that could exist

GREASE: Generally, every moving part should be greased. Our products leave the factory completely greased and lubricated, however, it is possible that with the time and usage of the product the elements degrease, either by loss of lubrication or by dirt. Clean and grease periodically the affected zones according to the manufacturer's specifications.

Verify that there are no loose, worn, or missing parts. Examine periodically every moving element to ensure that they are all tighten. Do not grease the platform guide because of its sliding dry system. **WARNING: Do not lubricate the support once connected to the electricity network.**

WEAR AREAS: Regularly examining the system's components to check if they are wearing out is a preventive measure that can reduce breakdowns. Verify possible lubricate leaks, grooves or bearings in poor condition.

MECHANICAL FIXING: Generally we call mechanical fixing to the components that are used to attach joint products, mainly screws and derivatives.

In some conditions of use, due to vibrations or impacts, many elements could lose their tightening torque or hold properties. Periodically examine that no elements are unleashed, especially in the moving parts of the set. Take into account and always respect the recommended tightening torque. Make sure to keep these areas free of water and moisture. Take special care not to get them wet when washing the set.

The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage.

Please note that cleaning must be carried out as described in this manual and that functionality must be verified before and after each use. Spencer Italia S.r.l. declines any responsibility for improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and compliance to Medical Device Directive 2017/745/EU.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with EU Regulation 2017/745.

12.3 PERIODIC OVERHAUL

The device must be overhauled every year by the Manufacturer, who uses internal and external technicians specialised and authorized by the Manufacturer itself.

In the absence of an above-mentioned overhaul, the device must be DECOMMISSIONED, as compliance with Regulation 2017/745/EU will cease and, despite the CE marking, the device will no longer meet the safety requirements guaranteed by the Manufacturer at the time of supply.

Spencer Italia srl declines any responsibility for incorrect operation or damage caused by use of devices that have not been regularly serviced.

The mattress and straps shall be replaced every two years.

Only overhaul activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

12.4 SPECIAL MAINTENANCE

Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l..

The end-user can replace only the spare parts indicated in § 15.

12.5 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase, which can be extended following annual overhauls.

Overhauls must be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer. In the absence of such annual overhauls, the device must be DISPOSED OF ACCORDING TO THE INFORMATION IN PARAGRAPH 15 AND NOTIFICATION MUST BE GIVEN TO THE MANUFACTURER.

The life span may be extended at the sole discretion of the Manufacturer or authorised centre if the safety requirements of the device are still met.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

13. TROUBLESHOOTING TABLE

COD.	ERROR'S DESCRIPTION	POSSIBLE CAUSE	POSSIBLE SOLUTION
1	VOLTAGE BATTERY IS LESS THAN 3V	- BATTERIES DISCONNECTED - BATTERIES DAMAGED	- BATTERIES PROPERLY CONNECTING - CHECK CABLE CONNECTION BATTERIES WITH ELECTRIC PLATE - REPLACE BATTERIES
2	VOLTAGE BATTERY IS LESS THAN 34V	OVERCHARGE DURING THE CHARGING BATTERIES OV	- COOL THE BATTERIES DURING 3/4 HOURS - USE THE CHARGER SUPPLIED BY SPENCER - REPLACE BATTERIES FOR NEW ONES
3	THE STRETCHER IS NOT RECEIVE VOLTAGE FROM AMBULANCE	- FAIL CONNECTION STRETCHER-SUPPORT - AMBULANCE'S BATTERY DISCONNECTED OR DAMAGED - POWER CABLES WRONG WAY CONNECTED BETWEEN SUPPORT AND AMBULANCE - ONE OF THE COPPER TERMINALS OF THE STRETCHER DID NOT CAME OUT CORRECTLY	- CHECK TERMINAL CONNECTION BETWEEN STRETCHER-SUPPORT - CHARGE THE BATTERIES OR REPLACE FOR NEW ONES - CHECK THE CABLE CONNECTION BETWEEN STRETCHER-SUPPORT-AMBULANCE - CHECK TERMINAL CONNECTION OF THE STRETCHER GOING INTO AND LEAVING THE ACCOMMODATION
4	VOLTAGE BATTERY IS LESS THAN 18V	BATTERIES ARE DISCONNECTED	CHARGE BATTERIES UNTIL ITS MAXIMUM LEVEL
5	THE STRETCHER GETS THE VOLTAGE THROUGH TERMINAL CONNECTION OF THE SUPPORT, BUT THEY HAVE NOT BEEN PRESSED	- THE STRETCHER IS NOT PROPERLY ENTERED IN THE SUPPORT - FAIL IN MICROSWITCH	- ENSURE THE STRETCHER HAS ENTERED WELL INTO THE SUPPORT AND THE CONNECTION TERMINAL ARE WELL CONENCTED - CHECK THE STATUS OF MICROSWITCH
6	BATTERIES ARE NOT DETECTED	- BATTERIES DISCONNECTED - BATTERIES CONNECTION DAMAGED - BATTERIES CHARGE VERY LOW	- ANCHOR WELL THE BATTERIES TO ITS BASE - REPLACE FOR NEW BATTERIES - CHARGE THEM WITH EXTERNAL SPENCER'S CHARGER
7	FAIL IN THE DRIVER OF ENGINE		CONTACT WITH TECHNICAL SERVICE
8	INTENSITY ABOVE 48V	- STRETCHER WITH OVERLOAD OF WEIGHT - BATTERIES UNCHARGED - THE STRETCHER NOT UP/DOWN GENTLY	- NOT LOAD THE STRETCHER WITH MORE THAN 300 KG - CHARGE THEM WITH EXTERNAL SPENCER'S CHARGER - SEARCH POSSIBLES SCRATCHS, INTERFERENCE OR LOCKED PIECES IN THE STRETCHER
9	VOLTAGE IS LESS THAN 15V	- STRETCHER WITH OVERLOAD OF WEIGHT - BATTERIES UNCHARGED - THE STRETCHER NOT UP/DOWN GENTLY	- NOT LOAD THE STRETCHER WITH MORE THAN 300 KG - CHARGE THEM WITH EXTERNAL SPENCER'S CHARGER - SEARCH POSSIBLES SCRATCHES, INTERFERENCE OR LOCKED PIECES IN THE STRETCHER
10	FAIL IN CHARGING PROCESS		CONTACT WITH TECHNICAL SERVICE

⚠ IMPORTANT CONSIDERATION

- The electrical system of the stretcher will automatically turn off in 15 minutes if there are not make any operation.
- In manoeuvre of up and down repeatedly to 20 cycles continuously, you should make a mandatory stop for at least 15 seconds for electro-hydraulic system of the stretcher recovers.
- It is recommended to remove the two batteries if the stretcher will spend long periods of inoperability.

- It is important to maintain a minimum of loading cycles every month to keep the good capacity of the batteries.
- Not leave without charge the batteries more than 2 MONTHS. It is possible they will stop working.

14. SPARE PARTS

If spare parts are required, please contact the manufacturer.

15. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

WARNINGS FOR THE CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/EU WEEE

At the end of his life, the product must not be disposed as household waste. Can be taken to special recycling centers provided by local government, or return it to the dealer on purchase of a new device of the same type and used for the same functions. Dispose of the product separately avoids possible negative consequences for the environment and human health resulting from inappropriate disposal and allows to recover the materials in order to obtain significant savings in energy and resources. The symbol on the label indicates separate collection of electrical and electronic equipment.

Warning: One incorrect disposal of electrical and electronic equipment could result in sanctions.

Notice

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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