AER0resc® - EASY Splint Universalschienen / Universal Splints

REF: HIMT02-SP-XL | HIMT02-SP-M | HIMT02-SP-XS | HIMT02-SP-XL-SET | HIMT02-SP-M-SET | HIMT02-SP-SET1 | HIMT02-SPT-XL | HIMT02-SPT-M | HIMT02-SPT-XS | HIMT02-SPT-XL-SET | HIMT02-SPT-M-SET | HIMT02-SPT-SET1

de Deutsch

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Achtung:

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

1. Safety Instructions

Attention:

- 1.1 The instruction for use does not replace detailed instructions / training of the user.
- 1.2 The AEROresc® - EASY Splint Universal Splints may only be used for the intended use described in these instructions for use.
- 1.3 Do not use for fractures of the spine and pelvis.
- 1.4 Do not use on open wounds. Open wounds must be covered with sterile dressings before using the splint.
- 1.5 Cutting the splint might result in sharp edges that need to be covered or folded.
- 1.6 For prolonged treatment duration (> 1 hour) a patient-side cushion pad should be used.
- 1.7 When fixing the splints, make sure that the blood flow is not interrupted and that the tips of fingers or toes are exposed to check the peripheral circulation.

4

- 5.2 The splints can be cleaned using mild cleaner or surface disinfectant. Please follow the instructions of the manufacturer.
- 5.3 Rinse with clear water to remove the cleaner and/or disinfectant.
- 5.4 Let the splints dry completely.

6. Technical Data

Dimensions and weight

Size XL: 100 x 11 cm | approx. 160 g

Size M: 50 x 11 cm | approx. 80 g

Size XS: 5 x 11 cm | approx. 4 g

Storage temperature: - 25 °C to 65 °C


Material

Coating: Polyethylene Foam (IXPE)

Core: Aluminum

7. Examples of use


Cervical Collar (Splinting)

Finger-Splinting

Tong Upper Arm-Splinting

Ankle-Splinting


Other applications are possible!


Double-layered Wrist-Splinting

Tong Forearm-Splinting

Knee- / Leg-Splinting

2. Intended Use

The AEROresc® - EASY Splint Universal Splints are used for immobilization and splinting of fractures and sports injuries (for example sprained ligaments). They can also be used as infusion splints.

3. Product description


3.1 The AEROresc® - EASY Splint Universal Splints are made of light and X-ray translucent materials. The coat is made of plastic foam (IXPE) - patient side „blue“ or „grey“ | Exterior „yellow“ or „olive green“ - and the core is made of an aluminum mixture.

3.2 The flexible material can be:

- cut to the required size with scissors (not serrated).
- preformed so that the affected body parts can be optimally stabilized.
- folded to increase stability, bent as required or combined with other splints.

4. Operation and handling

- 4.1 Form the splint so that a perfect fit is guaranteed to the body part.
- 4.2 Place the splint on the affected body part and secure it with elastic bandage or tape. For prolonged use, it is recommended to use a patient-side padding.
- 4.3 While fixing the splint ensure that you can control the blood circulation on fingers or toes.
- 4.4 Basic forms:

C-Curve

- ① To form the C-Curve, pull both side edges inwards.

Reverse C-Curve

- ② Form the C-Curve. Then form the side edges in the opposite direction.

T-Curve

- ③ Place the side edges on top of each other. Then bend half of each side until the bend has reached a 90 ° angle.


- 4.5 Examples of use see section "7. Examples of use".

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the relevant competent authority.

5. Cleaning and disinfection

- 5.1 Clean and disinfect the AEROresc® - EASY Splint Universal Splints after each use in accordance with the locally applicable hygienic guidelines.

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EC REP MedNet EC-REP GmbH
Borkstraße 10, 48163 Münster, Germany

Anping Guardian Medical Equipment Co., Ltd
500 meters, Nan Sucun, Anping County,
Hengshui City, Hebei Province, 053600, P.R.China

6

7