# READY HEAT

<b>EN</b>	<b>Ready-Heat™</b> Heated Disposable Temperature Management Blanket Instructions For Use .....	<b>2</b>
<b>DE</b>	<b>Ready-Heat™</b> Heizdecke zum Einmalgebrauch Gebrauchsanweisung .....	<b>3</b>
<b>FR</b>	<b>Ready-Heat™</b> Couverture de gestion de température jetable, chauffante Instructions d'utilisation .....	<b>4</b>
<b>IT</b>	<b>Ready-Heat™</b> Coperta autoriscaldante monouso per la regolazione della temperatura corporea Istruzioni per l'uso .....	<b>5</b>
<b>ES</b>	<b>Ready-Heat™</b> Manta de control de temperatura y calentamiento desechable Instrucciones de uso .....	<b>6</b>
<b>PT</b>	<b>Ready-Heat™</b> Cobertor de Aquecimento Descartável Instruções de Utilização .....	<b>7</b>
<b>SV</b>	<b>Ready-Heat™</b> Temperaturreglерande värmefilt för engångsbruk Bruksanvisning .....	<b>8</b>
<b>NL</b>	<b>Ready-Heat™</b> Verwarnde disposable temperatuurregelende deken Gebruiksaanwijzing .....	<b>9</b>
<b>DA</b>	<b>Ready-Heat™</b> Opvarmet engangstæppe til temperaturstyring Instruktioner for brug .....	<b>10</b>
<b>PL</b>	<b>Ready-Heat™</b> Podgrzewany jednorazowy koc termiczny Instrukcja użytkowania .....	<b>11</b>
<b>FI</b>	<b>Ready-Heat™</b> Lämmitetty kertakäyttöinen lämmönhallintapeite Käyttöohjeet .....	<b>12</b>
<b>EL</b>	<b>Ready-Heat™</b> Θερμαινόμενη κουβέρτα μίας χρήσης με ελεγχόμενες ρυθμίσεις θερμοκρασίας Οδηγίες χρήσης .....	<b>13</b>
<b>NO</b>	<b>Ready-Heat™</b> Oppvarmet temperaturkontrollteppe til engangsbruk Bruksanvisning .....	<b>14</b>
<b>CS</b>	<b>Ready-Heat™</b> Vyhřívaná přikrývka na jedno použití pro regulaci teploty Návod k použití .....	<b>15</b>
<b>TR</b>	<b>Ready-Heat™</b> Isıtmalı Tek Kullanımlık Sıcaklık Yönetim Battanitesi Kullanma Talimatı .....	<b>16</b>
<b>RU</b>	<b>Ready-Heat™</b> Одноразовое одеяло с подогревом и регулировкой температуры Инструкции по применению .....	<b>17</b>
<b>HU</b>	<b>Ready-Heat™</b> Egyszer használatos, fűtött testhőmérséklet-szabályozó takaró Használati utasítás .....	<b>18</b>
<b>KO</b>	<b>Ready-Heat™</b> 가온식 일회용 체온 관리 담요 사용 설명서 .....	<b>19</b>
<b>ZH<sub>hant</sub></b>	<b>Ready-Heat™</b> 可加熱一次性溫度可控毯子 使用說明 .....	<b>20</b>
<b>ZH<sub>hans</sub></b>	<b>Ready-Heat™</b> 可加热一次性温度可控毯子 使用说明 .....	<b>21</b>

## ENGLISH

---

### Ready-Heat™ Heated Disposable Temperature Management Blanket

This product is not sterile.

Made in the USA.

**Note:** This product is CE marked and manufacturer is ISO 13485 compliant.

**Description of Product:** Blue medical grade non-woven material with 9 oxygen activated heat packs.

**Indications and Usage:** To provide warm relief. Aids in the prevention of hypothermia.

**Contraindications:** See Warnings

#### Instructions for Use:

(Refer to Pages 22-23 for all Figures)

1. **DO NOT** open the outer protective package until ready to use, as activation will start immediately. Tear open the package at the **TEAR HERE** mark. *Figure 1.*
2. Remove the blanket from the outer protective package (keep outer protective package for blanket disposal). *Figure 2.*
3. Unfold the Ready-Heat Blanket to expose to air – heat activation will start automatically. *Figure 3.*
4. If the patient is not wearing a gown, place a drape over his/her body. **DO NOT PLACE THE BLANKET DIRECTLY ON EXPOSED SKIN.** *Figure 4.*
5. Place the blanket on the patient as indicated by the sticker, with the diagram sticker and the arrows pointing toward the head of the patient and facing upward. *Figure 5.*
6. Allow 15-20 minutes for the blanket to reach operating temperature, however the blanket will start warming immediately. *Figure 6.*
7. To access patient slide or turn the blanket. **DO NOT** fold blanket. *Figures 7 and 8.*
8. **FOR PROPER DISPOSAL:** After use fold the blanket, place back in the outer protective package and dispose with regular waste. If the blanket was in contact with blood-borne pathogens place back in the outer protective package and dispose with medical waste. *Figures 9 and 10.*

#### Warnings:

**\*\*Failure to follow these precautions may result in serious burns or other serious injuries.\*\***

1. Blanket is **NOT** sterile and is intended for single use only.
2. Careful supervision is required for use with the elderly, people with sensitive skin and those who can't feel heat or regulate heat. Do not use with infants and children.
3. Keep out of reach of children.
4. Monitor the temperature and cutaneous response of patients who are incapable of reacting, communicating and/or are without a sense of feeling every 10-20 minutes or according to institutional protocol.
5. Closely monitor vital signs for all patients on a regular basis.
6. Immediately notify physician of any vital sign instability.
7. **DO NOT** fold blanket while on the patient. Slide or remove blanket to access patient. Folding will cause heat retention magnification. *Figure 8.*
8. Remove blanket immediately if it becomes too hot or if the patient becomes too warm. Maximum temperature may reach 43°C.
9. **DO NOT** open the sealed pockets. **DO NOT** allow the contents of the sealed pockets to contact mouth, eyes or open wounds. If contact occurs, wash area with clean water. If contents are swallowed, call poison control center immediately.
10. **DO NOT** place directly on exposed skin.
11. **DO NOT** use under patient.
12. **DO NOT** use in MRI suite.
13. **DO NOT** use in hyperbaric oxygen chamber.
14. **DO NOT** put over operative field.
15. **DO NOT** place directly over transdermal medication patch.
16. **DO NOT** place instruments or objects on the blanket at any time. *Figure 11.*
17. **DO NOT** use on ischemic limbs.
18. **DO NOT** re-use
19. **DO NOT** allow blanket to come in contact with any plastic tracheal tubing.
20. **DO NOT** use blanket in the immediate area of any flammable anesthesia gases.

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TechTrade, LLC