
**Istruzioni d'uso**

**SPENCER 520**

IT

**Barella da campo**

**Instructions for use**

**SPENCER 520**

EN

**Field stretcher**

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**SPENCER 520**

DE

**Feldtrage**

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**SPENCER 520**

ES

**Camilla de campo**

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**SPENCER 520**

PT

**Maca de campo**

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|            |                                            |           |
|------------|--------------------------------------------|-----------|
| <b>1.</b>  | <b>MODELS .....</b>                        | <b>20</b> |
| <b>2.</b>  | <b>INTENDED USE .....</b>                  | <b>20</b> |
| <b>3.</b>  | <b>REFERENCE STANDARDS.....</b>            | <b>20</b> |
| <b>4.</b>  | <b>INTRODUCTION .....</b>                  | <b>21</b> |
| 4.1        | Using the manual .....                     | 21        |
| 4.2        | Device labelling and traceability .....    | 21        |
| 4.3        | Symbols.....                               | 22        |
| 4.4        | Warranty and service.....                  | 22        |
| <b>5.</b>  | <b>WARNINGS.....</b>                       | <b>22</b> |
| <b>6.</b>  | <b>SPECIFIC WARNINGS .....</b>             | <b>23</b> |
| 6.1        | Physical requirements for operators .....  | 24        |
| <b>7.</b>  | <b>RESIDUAL RISK .....</b>                 | <b>24</b> |
| <b>8.</b>  | <b>TECHNICAL DATA AND COMPONENTS .....</b> | <b>25</b> |
| <b>9.</b>  | <b>COMMISIONING.....</b>                   | <b>25</b> |
| <b>10.</b> | <b>OPERATING CHARACTERISTICS .....</b>     | <b>26</b> |
| <b>11.</b> | <b>PROPER USE .....</b>                    | <b>27</b> |
| <b>12.</b> | <b>CLEANING AND MAINTENANCE .....</b>      | <b>29</b> |
| 12.1       | Cleaning .....                             | 30        |
| 12.2       | Routine maintenance.....                   | 30        |
| 12.3       | Sheet replacement.....                     | 31        |
| 12.4       | Periodic overhaul.....                     | 31        |
| 12.5       | Special maintenance .....                  | 32        |
| 12.6       | Life span .....                            | 32        |
| <b>13.</b> | <b>TROUBLESHOOTING TABLE.....</b>          | <b>32</b> |
| <b>14.</b> | <b>ACCESSORIES.....</b>                    | <b>32</b> |
| <b>15.</b> | <b>SPARE PARTS.....</b>                    | <b>32</b> |
| <b>16.</b> | <b>DISPOSAL .....</b>                      | <b>32</b> |

## 1. MODELS

The following basic model may be subject to implementation or change without notice.

- SPENCER 520

## 2. INTENDED USE

Field stretchers are devices for lifting and transporting people, also suitable for use in field rescue situations. If they do not have supports, they can be stored using the special Porta Stock stand. The patient is not expected to intervene on the devices.

### Target patients

There are no particular indications related to the patient group. The product configuration is able to accommodate any subject as long as he/she is within the maximum capacity of the device.

### Patient selection criteria

The patients foreseen are those with injuries that prevent them from walking in a given rescue situation, or who are in a state of unconsciousness. No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

### Contraindications and unwanted side effects

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

### Users and installers

The intended users are field rescue workers, typically nurses, stretcher bearers and doctors present in the field.

These devices are not intended for use by lay people.

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands.

**Operators must be able to provide the necessary patient care.**

**Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device.**

**For patient loading techniques of particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (in addition to 2 as required in standard conditions).**

The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.

### User training

*Note: despite all efforts, laboratory tests, trials, and instructions for use, standards do not always reproduce practice, so the results obtained under actual conditions of product use in the natural environment could sometimes differ significantly.*

*The best instructions are the continuous practice of use under the supervision of competent and trained personnel.*

- Regardless of your level of past experience with similar devices, the user should carefully read and understand the contents of this manual before installing, operating, or servicing this product. In case of any questions, please contact Spencer Italia S.r.l. for the necessary clarifications.  
The product must be used only by personnel trained in the use of this specific product. The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. **This documentation must be kept for at least 10 years after the end of the product's life and must be made available to the competent authorities and/or the Manufacturer when requested. In the absence of such documentation, the relevant bodies will apply any foreseen sanctions.**
- Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.
- The product must be put into use only by personnel trained in the use of this product and not on other similar products.

*Note: Spencer Italia S.r.l. is always at your disposal to organise product training.*

## 3. REFERENCE STANDARDS

As a Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.

| REFERENCE              | TITLE OF DOCUMENT                                                                                       |
|------------------------|---------------------------------------------------------------------------------------------------------|
| EU REGULATION 2017/745 | EU REGULATION 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices |

## 4. INTRODUCTION

### 4.1 Using the manual

The purpose of this manual is to provide healthcare professionals with the information necessary for safe and appropriate use and maintenance of the device.


*Note: The Manual is an integral part of the device and therefore it must be kept for the entire life of the device and must accompany it in any changes of use or ownership. If any instructions for use for products other than the one received are present, please contact the Manufacturer immediately before use.*

The instructions for use of Spencer products can be viewed and downloaded from [www.spencer.it/risorse-media/manuali-uso](http://www.spencer.it/risorse-media/manuali-uso) on the Spencer website or you can contact the Manufacturer. Exceptions are those items whose essentiality and reasonable and predictable use are such that it is not necessary to draw up instructions, in addition to the following warnings and indications on the label.

Regardless of your level of experience with similar devices in the past, it is advisable to carefully read and understand the contents of this manual before installing, operating, or servicing this product.

### 4.2 Device labelling and traceability

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

In the event of damage or removal, request a duplicate from the Manufacturer, or else the warranty will be void as the device can no longer be traced.

Regulation 2017/745/EU requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly by Spencer Italia S.r.l., inform After Sales (see § 4.4).

### 4.3 Symbols

| Symbol                                                                            | Meaning                                          |                                                                                                                 |                                                                                                                      |
|-----------------------------------------------------------------------------------|--------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
|  | Device in compliance with EU Regulation 2017/745 |                                | See the user manual.                                                                                                 |
|  | Medical device                                   |                                | Serial Number                                                                                                        |
|  | Manufacturer                                     |                                | Product code                                                                                                         |
|  | Date of manufacture                              |                                | Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death. |
|  | Unique Device Identifier                         | Production identification<br>Alphanumeric code that identifies the production units of the device, composed of: |                                                                                                                      |
|  |                                                  |                                                                                                                 |                                                                                                                      |
|                                                                                   |                                                  | (01)0805771123XXXX                                                                                              | Company prefix (XXXX=GTIN)                                                                                           |
|                                                                                   |                                                  | (11) YYMMDD                                                                                                     | Date of manufacture                                                                                                  |
|                                                                                   |                                                  | (NN) 1234567890                                                                                                 | (NN)1234567890 NN=10 =>LOT/<br>NN=21=>SN                                                                             |

### 4.4 Warranty and service

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

For information regarding correct interpretation of the instructions for use, maintenance, installation or return, please contact Spencer Customer Service tel. +39 0521 541111, fax +39 0521 541222, e-mail [service@spencer.it](mailto:service@spencer.it).

To facilitate service, always indicate the lot number (LOT) or serial number (SN) on the label attached to the package or device itself.

**Warranty and service conditions are available at [General Terms and Conditions | Spencer](#).**

*Note: Record and keep with these instructions: lot (LOT) or serial number (SN), if present, place and date of purchase, date of first use, date of checks, user name and comments.*

## 5. WARNINGS

Warnings, notes, and other important safety information are provided in this section and are clearly visible throughout the manual.

At least every 6 months, it is important to check for updated instructions and any changes involving your product. This information is freely available on the website [www.spencer.it](http://www.spencer.it) on the specific product page.

### Product features

**Use of the product for any purpose other than that described in the User Manual is prohibited.**

- Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device are detected, the device must be immediately removed from service and the manufacturer must be contacted.
- The use by untrained personnel may result in injury to the patient, rescuer or others.
- Inadequate disinfection procedures may involve cross-infection risks.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruptions.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- During use, position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids, when applicable.
- For devices intended for the transport of patients, always respect the maximum load capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total loading weight of the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
- For devices intended for the transport of patients, make sure, before lifting, that operators have appropriate physical conditions, as

indicated in the manual.

- The maximum weight supported by each rescuer must comply with the requirements prescribed by the law of the country of use, in the field of Health and Safety at Work.
- **The warranty seals, where present, must not be removed;** in such case, the manufacturer will no longer recognize the product warranty and will accept no responsibility in case of incorrect operation or damage caused by the product itself.
- Avoid contact with sharp objects.
- Operating temperature: from -10°C to + 50°C

#### Storage

- The device should not be exposed to or come into contact with any source of combustion or inflammable agents. Store in a cool, dry, dark place and do not expose to direct sun.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store and transport the device in its original packaging. Failure to do so makes the warranty void.
- Storage temperature: from -20°C to +60°C

#### Regulatory requirements

**As Distributor or End-User of the products manufactured and/or marketed by Spencer Italia S.r.l., users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory.**

- Promptly and in detail notify Spencer Italia S.r.l. (immediately in the quotation request phase) about possible fulfilments by the Manufacturer necessary for the compliance of products with specific legal requirements of the territory (including those deriving from regulations and/or regulatory provisions of another nature).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, exactly as specified in the relevant user's manual.
- **Actively contribute to safety checks on product** sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- The distributor or final user is aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damage that might occur. Therefore Spencer Italia S.r.l. expressly disclaims any responsibility and/or liability for your non-compliance with the present regulatory provisions.
- With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the Ministry of Health, within the terms and in the manner established by one or more ministerial decrees, notifying also the Manufacturer. Public or private health care professionals are required to notify the Manufacturer of any other incident that may allow the adoption of measures to ensure the protection and health of patients and users.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with Regulation 2017/745/EU Medical Devices.

#### 6. SPECIFIC WARNINGS

- Follow the procedures approved by the Emergency Medical Service for the immobilization and transport of patients.
- Follow the procedures approved by the Emergency Medical Service for the positioning and transport of patients.
- Do not use if the device or its parts are pierced, torn, frayed or excessively worn out.
- Make sure, before lifting, that the operators have a firm grip on the device.
- When lifting the device in stretcher mode, grab the telescopic handles at the head side and the ones at foot side which are placed on the main frame. **Don't use the handles placed on the seat frame to lift the device.**
- Avoid pulling the device on rough surfaces.
- Do not lift by crane or other mechanical lifters.
- Do not use drying machines.
- The device is intended for recovery and transfer of patients and the patient should not be left permanently on the device
- Do practice with an empty stretcher/chair, to make sure you become familiar with the manoeuvres.
- For patient loading techniques, for particularly heavy patients, for working on uneven ground or in special and unusual circumstances, the presence of more operators is recommended (not only 2 as required under standard conditions).
- Before each use, check the integrity of the belts and their hooks, as specified in the user's manual. In case of malfunction or damage that may compromise the functioning and safety of the device, patient or operator, it is necessary to replace the belts.
- Make sure the belts are properly fastened to the frame of the stretcher.
- Always immobilize the patient, using the straps supplied by the manufacturer; lack of immobilization could cause serious damage.
- To optimize the immobilization of the patient, in addition to the belts provided, the use of a thorax belt is required.
- Make sure the sheet is properly secured to the frame of the stretcher.

- Make sure that no sheet interferes with the folding system, with the handle movement and with any other operation on the device.
- Do not operate if the weight has not been distributed correctly.
- Always grasp only the handles to carry the stretcher/chair and not other parts of the device.
- Avoid excessive force when changing the configuration of the device. Excessive force can compromise the integrity of the product.
- Before folding the device for storage, check that it is fully extended. The partial extraction of rear handles, could cause damage to the device.
- Don't try to carry out configuration changes if patient is on the device.
- Before configuration change, make sure that nobody has hands placed between sliding or moving systems. It could lead to injury.
- Hold the stretcher firmly still when the patient is on it.
- Partial extraction of telescopic handles can cause its sudden closing or opening, with related risks for patient and operators. Always make sure that the handles are fully inserted.
- Failure to use safety belts can cause serious risks for the patient. Always use the belts and make sure that patient is properly secured.
- A sheet not properly fixed, can detach on loading causing risk to the patient. Always ensure a proper installation of the sheet.
- Contact with sharp objects or abrasive surfaces, may compromise the integrity of sheet and belts, causing risks for the patient.
- Failure to verify the integrity of belts and sheets, can result in breakage compromising patient safety.
- Failure to comply with warnings for operators, may cause danger of crushing cause by handling mechanisms.
- When opening the device from folded position used for storage, ensure that nobody has hands placed near the folding point; crushing injuries may occur.
- Pay a lot of attention to possible obstacles (water, ice, debris, etc.) on the route of the stretcher/chair, because they could cause loss of balance of the operator and compromise the proper functioning of the device. If the path ~~free~~ cannot be made free of obstacles, choose an alternative path.
- Water, ice or debris on the stairs can cause loss of balance for the operator and jeopardize the proper functioning of the device. In order to avoid injury, clear the way of obstacles or choose an alternative route.
- Do not alter or modify the stretcher arbitrarily to adapt it to your own purposes. The modification may cause unforeseeable functioning and damage to the patient or operators. In any case the warranty will be void, relieving the manufacturer from any liability.
- Always verify that the configuration change has been carried out completely by reaching the respective end of stroke.
- It is preferable not to use the device in chair mode in the event of suspected cervical or spinal trauma, or fractures.
- For a safe use of the device, the presence of at least two operators is necessary.
- Don't use the wheels to go up or down stairs. If there are obstacles, the device must be lifted.

## 6.1 Physical requirements for operators

The device is intended for professional use only. All operators must be trained to transport patients safely and efficiently. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others.

Operators using the device should be physically able to use the device and have good muscle coordination, as well as strong back, arms, and legs for lifting and holding and be able to grasp the device firmly with both hands.

**Operators must be able to provide the necessary patient care.**

**Users must be able to safely lift and handle the weight of the stretcher and patient together and any other equipment used with the device.**

**For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of multiple operators is recommended (not just 2 as expected in standard conditions).**


The ability of each operator must be assessed before the definition of rescuer roles in use of the stretcher.

## 7. RESIDUAL RISK

No residual risks have been identified, i.e. risks that could arise despite compliance with all the warnings contained in these instructions for use.

## 8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without prior notice.

| N° | Description        | Material     | N° | Description                   | Material       |
|----|--------------------|--------------|----|-------------------------------|----------------|
| 1  | Telescopic handles | Steel/PVC    | 7  | Stretcher's rear handles      | PVC            |
| 2  | Backrest sheet     | PVC          | 8  | Configuration change crossbar | Steel          |
| 3  | Seat sheet         | PVC          | 9  | Seat frame                    | Steel          |
| 4  | Seat frame handles | PVC          | 10 | Main frame                    | Aluminum/Steel |
| 5  | Footrest           | Steel        | 11 | Chair's rear handles          | Steel/PVC      |
| 6  | Wheels             | Polyurethane | 12 | Ground support feet           | Steel          |

| CHAIR MODE                        | Size                |                     | Stretcher mode                        | Size | GENERAL                            |            |
|-----------------------------------|---------------------|---------------------|---------------------------------------|------|------------------------------------|------------|
| Total height                      | 1305                |                     | Length with telescopic handles closed | 1945 | Width (mm)                         | 545        |
| Seat/Backrest height              | 480/825             |                     | Length with telescopic handles opened | 2055 | Weight (kg)                        | 14,2 ± 0,5 |
| Seat width                        | 380                 |                     |                                       |      | Loading capacity <sup>1</sup> (kg) | 160        |
| Depth                             | Rear handles closed | Rear handles opened | Lying part height                     | 100  |                                    |            |
|                                   | 805                 | 945                 |                                       |      |                                    |            |
| Dimensions folded 1100x545x170 mm |                     |                     |                                       |      |                                    |            |

<sup>1</sup> Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total load of the weight on the product, the operator must consider the weight of the patient, the equipment and the accessories.

NOTE: Dimensions are subjected to tolerances of ± 10mm

## 9. COMMISSIONING

For first use, check that:


- Packaging is intact and has protected the device during transportation
- Check that all parts included in the packing list are present
- General functionality of the device
- If present, check that the patient bed adjustment systems lock in place properly
- If present, check that all controls are working properly
- If present, check that all unlocking mechanisms are working properly

Do not modify the stretcher structural, leverage or traction parts in any way as this could cause damage to the patient and/or the rescuers. Failure to take the above measures will preclude safe use of the device, resulting in risk of damage to the patient, operators and the device itself.

**For subsequent use, perform the operations specified in paragraph 12.**

If the above conditions are met, the device may be considered ready for use; otherwise, you must immediately remove the device from service and contact the Manufacturer.

## 10. OPERATING CHARACTERISTICS

| Element | Description                   | Function                                                                                                                                           |
|---------|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| A       | Telescopic handles            | They can be used both in chair and stretcher mode, for transport and lifting of the device. They can be extracted to make these operations easier. |
| B       | Seat and backrest sheets      | They are the patient support. They are fixed to the frame by means of belts with quick-release buckles.                                            |
| C       | Seat frame handles            | Used to lift the device when in chair mode, by the operator positioned at the front side                                                           |
| D       | Wheels                        | They make small movements or manoeuvres of the device easier when used in chair mode                                                               |
| E       | Footrest                      | Allows the patient to rest the feet when the device is in chair mode and restricts longitudinal movements when used in stretcher mode              |
| F       | Stretcher rear handles        | They are the terminal component of the movable part of the main frame. They are used to lift the device in stretcher mode                          |
| G       | Configuration change crossbar | Used to carry out configuration changes                                                                                                            |
| H       | Seat frame                    | It supports the patient's weight when the device is in chair mode. It's the moving component that allows the configuration change                  |
| I       | Main frame                    | The main structure of the device                                                                                                                   |
| L       | Chair's rear handles          | Used by the operator placed at the rear side to lift the device when in chair mode                                                                 |
| M       | Ground support feet           | They support the device in stretcher mode when placed on the ground                                                                                |

## 11. PROPER USE

Before transferring, lifting or transporting the patient, primary medical evaluations have to be performed. Once the diagnosis has been assured, it is preferable (if possible) to ask the patient to actively collaborate during transfer onto the stretcher making sure the patient is fully aware of all risks. Before loading the patient, place the device near the patient; always maintain the safety belts fastened during the manoeuvres.

### Belts positioning

Both in chair and stretcher mode, is essential to use the supplied safety belts. The device comes with 3 belts.

The longer belt, intended to immobilize chest, should be positioned immediately above the ground support feet, passing it behind the backrest as shown in thumbnail **A**.


The location of this belt can be different if used with particularly small patients, for which it could be necessary to place the belt at a lower height. **Never pass this belt around the neck!**

The seat belt, must be passed under the seat frame.

The shorter belt, intended for leg immobilization, must pass between the vertical support frame of the seat and the seat as shown in the thumbnail **B**. When the device is used in stretcher mode, the belts should be placed in the same position.

**In order to allow faster rescue operations, is recommended to maintain the belts installed on the device and to open them only when loading the patient on the device.**

After patient positioning, fasten the belts and tighten them to ensure a correct immobilization.

### Telescopic handles

The telescopic handles have been made to facilitate lifting of the device when in stretcher mode or to make tilting easier when in chair mode. They are kept in position by a mechanism that is inserted automatically when reaching the excursion limits both during opening and closing.


To unlock the mechanism allowing handle movement, press downwards and move simultaneously.

Pull the handles until they are fully open. It is essential to verify that the locking system has been properly inserted.

To do this, repeatedly try to slide the handles without applying pressure downwards. **There must be no movement.**


**DO NOT LIFT THE DEVICE IF THE HANDLES ARE NOT COMPLETELY EXTRACTED AND LOCKED.**

Similarly it is possible to close the telescopic handles verifying that they are locked in closed position.

### Lifting the stretcher with patient

- Secure the patient to the stretcher with the supplied belts adjusting properly the tension according to the patient's clinical status. The operators should be at least two, one at the head side and the other one at the foot side.
- Using the appropriate lifting technique to minimize the effort for the operator, grab the handles on each side and raise the stretcher keeping it horizontal.
- For easier lifting, is possible to extract the handles as described in the previous paragraph.
- For lifting operations is strictly forbidden to grasp the handles of the seat frame. Lifting the device using such handles could cause its sudden folding.**

## Configuration change


During changes of configuration, the patient must not be on the device. It is recommended a preliminar evaluation intended to choose the best configuration.

### Stretcher Mode -> Chair Mode

The configuration change between stretcher mode and chair mode, takes place in these two steps:


1. Raise the device from the telescopic handle's side, and bring it up until vertical position is reached.
2. Keeping the device vertical, place a foot under the configuration change crossbar to stabilize the device and press down until the main frame has reached the end of stroke. If the operation has been carried out successfully, the distance between the top of the handle and the start of the main frame will be 2,5 cm.

### Chair Mode -> Stretcher Mode

This configuration change is performed similarly to the previous one:

1. Place the device in chair mode, bringing the footrest on the ground.
2. Place a foot on the configuration change crossbar.
3. Lift up the device by the telescopic handles or the main frame.
4. When the device is fully extended, bring it to the ground.

**When changing configuration, make sure that nobody has placed the hands near sliding systems or moving parts because it could result in crush injuries.**

## Transport in Chair Mode

The device can be used in chair mode with the wheels on the ground for manoeuvres in restrained space in order to reach a position suitable for lifting and transport.


Use the device, with or without the patient on it, using only the rear wheels, is not permitted.

The transport by lifting in chair mode, is particularly indicated in confined spaces such as stairs.


For this operation the presence of more operators is required, at least one at the front side and at least one at the rear side.

Proceed as follows:

- Ensure that patient's clinical conditions are suitable for transport in sitting position.
- Place the device in chair mode as described above verifying that the configuration change has been carried out properly.
- Open the rear handles until they are fully opened.
- **Place the patient on the device and secure it with provided belts.**
- **Coordinating movements, lift the device using the rear handles and the handles of the seat frame.**
- **Handling on wheels, is useful for short moves or to facilitate manoeuvres prior to subsequent lifting.**
- **Ensure that patient's clinical conditions are suitable for transport in sitting position.**

- Place the device in chair mode as described above verifying that the configuration change has been carried out properly.
- Grasp the telescopic handles from the rear side, place a foot on the configuration change crossbar in order to limit movements of the chair and tilt the device towards you.
- In order to limit risks for the patient, especially if heavy, this operation should not be performed by a single operator.
- **Move the device sliding it on the ground using the wheels.**
- **For obstacles higher than 1cm, is necessary to lift the device.**

#### Device folding




In order to minimize storage space, the device can be folded as follows:

1. Bring the device in **Stretcher mode**, making sure that it is fully extended
2. Lay it on a side end make sure that rear handles are closed.
3. Move to the back of the stretcher and fold the two halves towards you until they are in contact. If this operation is difficult, check that the seat frame is fully extended, allowing the full extension of the rear handles of the stretcher.  
To keep the stretcher in closed position preventing unwanted opening, the supplied belts can be used, wrapping them around the folded device. Adjust the belt's tension so that the two halves of the stretcher are in contact and stable.

4. To return the stretcher to extended position, carry out the same operations in reverse order.

**Warning, crushing risk! Don't place hands near the folding point.**

## 12. CLEANING AND MAINTENANCE

As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the final destination Country (including laws and norms regarding technical specifications and/or safety requirements) of the goods and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.

**Spencer Italia S.r.l. declines all responsibility for any direct or indirect damage which is the consequence of improper use of the product and spare parts and/or in any case of any repair carried out by a person other than the Manufacturer, who uses internal and external technicians authorised; to do so; moreover, doing so will invalidate the warranty.**

- The operator must always wear adequate personal protection such as gloves and mask etc. during all checking, maintenance and cleaning procedures.
- Establish a maintenance program and periodic testing, identifying an employee responsible for overseeing. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in the user's manual.
- The frequency of inspection is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.
- Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, that is the result of incorrect repairs or use of products made by Spencer Italia S.r.l. Repairs must necessarily be carried out by an authorized Spencer Italia service centre, which in using original spare parts will provide a quality repair service in strict accordance with the technical specifications given by the manufacturer. Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, which is a result of improper use of spare parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to repair or make substitutions on this product and parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to do so; the warranty will also be invalidated.

- Use only original components, spare parts and or accessories, approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device.
- For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.
- All maintenance and revision must be recorded and documented with the corresponding report for technical assistance; documentation shall be maintained for at least 10 years from the end of life and must be made available to the competent authorities and/or the manufacturer if requested.
- The cleaning schedule for reusable products must be performed in accordance with the directions provided by the manufacturer in the user manual, in order to avoid the risk of cross-infection due to the presence of secretions and/or residuals.
- The device and all its components, after washing, should be allowed to dry completely before storing.
- If required, lubrication must be carried out after cleaning and complete drying.

## 12.1 Cleaning

**Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals. The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.**

The exposed metal parts are usually treated and/or painted in order to increase their resistance to external agents. Clean the exposed parts with water and delicate soap; **never use solvents or stain removers.**

Rinse thoroughly with warm water making sure that you have removed all traces of detergent, which could degrade or compromise the integrity and durability of the device. **The use of high pressure water should be avoided** because water penetrates the joints creating the risk of corrosion of components. Allow the device to dry thoroughly before storing. Drying after washing or after use in wet environments must be natural and not forced; do not use flames or other sources of direct heat.

If you need to carry out **disinfection**, use products that do not have corrosive or solvent action on the materials constituting the device. Be sure to take every precaution to ensure that there is no risk of cross-infection or contamination for patients and / or operators.

To maintain the polished appearance of the frame parts we recommend the use of polish cleaner Spencer STX 99 or in alternative creams or waxes used for polishing car bodywork.

## 12.2 Routine maintenance

**Establish a maintenance program and periodic testing routine, identifying an employee responsible for this. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in following paragraphs are inspected.**

**All maintenance and periodic servicing activities must be registered and kept together with the servicing reports. These documents have to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.**

Routine maintenance of the device must be carried out by operators in possession of specific qualifications, trained and experienced in the use and maintenance of the device.

**The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.**

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Correct fixation of all nuts, bolts and screws
- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- The seat or the backrest has no structural damage or cracks
- Welded areas are intact, without any cracks or breaks
- The seat belts, sheets, moving parts, wheels and handles are intact and functioning properly.
- Lubrication of moving parts
- State of use of wheels and breaking system
- The wheels are correctly fixed, they are stable and turn properly
- The wheels are free from debris.
- The device opens and locks properly
- The device opens and closes properly
- There are belts for the immobilization of the patient and they are intact and functioning.
- Posterior telescopic handles can be positioned, folded and locked properly in all positions as described in this manual.
- Handles are stable and have no risk of slipping out
- The locking mechanism of the telescopic handles works properly

**The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.**


The only part to be lubricated if necessary and at least once a month, is the rotation axis of the wheels. For this maintenance, no disassembling is needed.

Lay the device on a side, put a drop of lubricant in the inner part like shown in the picture and turn the wheel repeatedly in order to distribute the lubricant over the entire inner surface.

Do the same for the other wheel.

Before lubricating, be sure to remove any dust or dirt in the areas, if necessary use compressed air.

Please note that cleaning must be carried out as described in this manual and that functionality must be verified before and after each use. Spencer Italia S.r.l. declines any responsibility for improper functioning or damages caused to the patient or user by the use of devices not subject to routine maintenance and will void the warranty and compliance to Medical Device Regulation 2017/745/EU.

**!** Every 2 years starting from the purchase date, is necessary to replace sheets and patient belts.

### 12.3 Sheet replacement



Sheet replacement can be carried out easily, thanks to the quick release buckles.

Place the device in stretcher mode and bring it to the ground with the quick release buckles facing upward. Remove old sheets releasing all the buckles.

Install the new sheet making them pass around the red areas identified by the numbers like shown in the picture. In the parts **1** and **2** the upper and lower side of the backrest sheet must be fixed (**A**).

In the parts **3** and **4** respectively the widest part and tighter one of the seat sheet must be fixed (**B**)

Do not pass the sheets around other part of the frame because movement of the device will be compromised. After the correct installation of the sheet has been verified, adjust the strap to ensure that the sheets are tight.

Use only original or Spencer Italia S.r.l. approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Otherwise, we decline all responsibility regarding incorrect operation or any damage caused by the device to the patient or the operator, invalidating the warranty and invalidating compliance with Regulation 2017/745/EU Medical Devices.

### 12.4 Periodic overhaul

**The device must be overhauled every year by the Manufacturer, who uses internal and external technicians specialised and authorized by the Manufacturer itself.**

**In the absence of an above-mentioned overhaul, the device must be DECOMMISSIONED**, as compliance with Regulation 2017/745/EU will cease and, despite the CE marking, the device will no longer meet the safety requirements guaranteed by the Manufacturer at the time of supply.

Spencer Italia srl declines any responsibility for incorrect operation or damage caused by use of devices that have not been regularly serviced.

The mattress and straps shall be replaced every two years.

The gas spring on the front legs must be replaced every 2 years. More heavy-duty or particular conditions of use may, however, require it to be replaced earlier. This need is evaluated by the manufacturer or by the service centre that carries out the overhaul.

Failure to replace it within the prescribed period could compromise the safety of the stretcher.

Only overhaul activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.

## 12.5 Special maintenance

**Special maintenance can only be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer itself.**

Only maintenance activities carried out by specialised technicians authorised by the Manufacturer are considered valid by Spencer Italia S.r.l.. The end-user can replace only the spare parts indicated in § 15.

## 12.6 Life span

**The device, if used as described in the following instructions, has a life span of 5 years from the date of purchase, which can be extended following annual overhauls.**

Overhauls must be carried out by the Manufacturer, who uses internal and external technicians specialised and authorised by the Manufacturer. **In the absence of such annual overhauls, the device must be DISPOSED OF ACCORDING TO THE INFORMATION IN PARAGRAPH 16 AND NOTIFICATION MUST BE GIVEN TO THE MANUFACTURER.**

The life span may be extended at the sole discretion of the Manufacturer or authorised centre if the safety requirements of the device are still met.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

## 13. TROUBLESHOOTING TABLE

| PROBLEM                                                      | CAUSE                                                                                 | SOLUTION                                                                                                                                                        |
|--------------------------------------------------------------|---------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| It is not possible to carry out the configuration            | Sliding systems or other moving parts are damaged                                     | Put the device immediately out of service and contact the service centre                                                                                        |
|                                                              | After replacement of the sheets, they have been fixed around wrong areas of the frame | Verify that the new sheets are installed properly                                                                                                               |
| Wheels do not allow an easy sliding                          | Dirt, debris and absence of lubrication limit the rotation of the wheels              | Carry out cleaning and lubrication procedures described in paragraph 12.2. If the problem persists, put the device out of service and contact the manufacturer. |
| The sheet doesn't support the patient's weight               | The quick release buckles are not properly engaged                                    | Verify the proper closure of the buckles placed on the rear side of the sheet                                                                                   |
|                                                              | Worn cloth                                                                            | Replace the sheet as described in the manual                                                                                                                    |
| Sheets are not tight enough                                  | The straps are loose due to use                                                       | Verify that the belts are tight enough                                                                                                                          |
|                                                              | Worn cloth                                                                            | Replace the sheet as described in the manual                                                                                                                    |
| The frame is deformed                                        | The device has been subjected to excessive load or has been used improperly           | Put the device immediately out of service and contact service center                                                                                            |
| The telescopic handles don't lock in closed or open position | The telescopic handles have not been completely extracted or inserted.                | Extract or insert the telescopic handles until the end of the stroke is reached                                                                                 |
|                                                              | The locking system is damaged.                                                        | Put the device immediately out of service and contact the service center                                                                                        |
| The rear handles are unstable or deform under loading        | Handle assembly damaged                                                               | Put the device immediately out of service and contact the service centre                                                                                        |

## 14. ACCESSORIES

There are no accessories for this device

## 15. SPARE PARTS

|         |                                  |
|---------|----------------------------------|
| RIST001 | BLACK WHEEL Ø 150mm WITH BEARING |
| RIST019 | SHEETS KIT SPENCER 520           |
| RIST020 | BELTS KIT SPENCER 520            |

## 16. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

---

**Notice**

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change.

The images are included as examples and may vary slightly from the actual device.

**© Copyright Spencer Italia S.r.l.**

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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