# RUSSELL CHEST SEAL®

**Intended Use:**

The Russell Chest Seal® is designed for the immediate management of an open chest wound and to minimise the risk of developing tension pneumothorax

**1 emergency chest seal w/integrated valve.  
20cm x 15cm (8" x 6")**

Medical Device

Sterilized using irradiation

Do not re-use

Do not use if package is damaged

Does not contain a presence of natural rubber latex

Do not re-sterilize

Consult instructions for use

Importer

Use by date

Single sterile barrier system

Safeguard Technologies Ltd.  
Willow Grove, Delgany,  
Co Wicklow, A63 XY89, Ireland  
Authorized Representative in the EC

Unique Device Identifier

Catalogue number

Batch Code

Safeguard Medical Technologies Ltd.  
The Old Rectory,  
Hope Under Dinmore,  
Herefordshire, HR6 0PW  
United Kingdom

[www.safeguardmedical.com](http://www.safeguardmedical.com)

NSN:6510-99-490-5137

**Intended use:**

The Russell Chest Seal® is designed for the immediate management of an open chest wound and to minimise the risk of developing tension pneumothorax.

**User profile:**

The Russell Chest Seal® is designed to be used by medical professionals in the pre-hospital setting.

**Patient profile:**

The Russell Chest Seal® is designed to be used on those patients with an open chest wound.

**Clinical benefits:**

- Valve mechanism allows air to escape from the pleural cavity without being drawn back in.
- Minimises risk of developing tension pneumothorax and lung collapse.

**Performance characteristics:**

- Fast, easy and intuitive to apply.
- Patented hydrogel base.
- 4 large peripheral apertures allow drainage of air and blood.
- Large central hole in base material allows quick accurate placement with reduced risk of obstructing the wound opening.
- Low profile and highly flexible.
- Compact: easy to pack into medical kit.


**Product claims:**

- The device is comfortable for the patient during use.
- The device will remain adhered over many hours.
- Latex free.

**Warnings and residual risks:**


- Do not place over protruding object as this may worsen injury or damage valve.
- Excessively wet or greasy skin will reduce the ability of the adhesive to function. If patient's chest is suspected to be saturated with lipid-based lotion, cleanse with a skin cleansing wet wipe prior to application.
- Do not smother the valve as this will risk loss of functionality.
- The device has been designed for short term use only, until such time as proper medical/surgical intervention can be initiated. For use up to 24 hours.
- Do not cut the dressing as this will risk loss of functionality.
- In case of intolerance, allergic reaction, or irritation to the dressing, remove and clean the area carefully, and seek medical advice.
- Device is single-use only. Reuse will risk cross-infection and loss of product functionality.
- Use past shelf life risks cross-contamination or loss of functionality.
- Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use if sterile barrier is compromised.
- Dispose of immediately after use, as clinical waste as per local protocols. Improper disposal risks cross-contamination.
- If difficulty is experienced on removing the dressing, it should be irrigated or soaked with water or sterile saline solution.
- Do not obstruct the valve, as this may reduce or stop its ability to function properly.
- Any serious incident that has occurred in relation to this device should be reported to **manufacturer (vigilance@safeguardmedical.com) and applicable national competent authority.**

**1**

Remove any loose debris or fluid from the skin around the wound. Note that excessively wet or greasy skin will reduce the ability of the adhesive to function.


**2**

Remove the dressing from its packaging and unfold the dressing. Line up the dressing so that the transparent central portion of the valve is positioned directly over the wound.

Remove one of the plastic butterfly liners and place the exposed adhesive surface onto the patient.

**3**

Remove the second liner, and gently firm down across the dressing's adhesive area to ensure that the dressing is properly attached.

DO NOT obstruct the valve, as this may reduce or stop its ability to function properly. Frequently monitor the valve to ensure that it is properly functioning.

## Management of open pneumothorax and role of the Russell Chest Seal®

Treatment of open pneumothorax involves applying a seal to the chest which – by means of a valve – will allow blood and air to be expelled (during exhalation) and prevent air from re-entering the chest on inhalation, allowing a negative pressure to be created and the lung to re-inflate.

The Russell Chest Seal® is a hydrogel-based dressing with an integral low-profile, four apertures, valve. The base material of the chest seal is a hydrogel matrix, which comprises a water-containing gel, mounted onto a non-woven fabric with silicon coated release liners.

The Russell Chest Seal® is used to cover penetrating chest wounds, and the valve mechanism allows air to escape from the pleural cavity without being drawn back in. This aims to minimise lung collapse associated with this injury and to facilitate lung re-inflation. The device will remain adhered over many hours and is comfortable for the patient in prolonged use.

On inhalation, even a small negative pressure will make the valve close, even in the presence of fluid such as blood. During repeated respiratory cycles, air can pass out through the seal with each exhalation and the lung will expand a little more with each inhalation. Gradually, the lung will re-inflate and oxygenation should improve.

The Russell Chest Seal® has no latex parts and is designed to remain highly occlusive in order to make an effective seal. By incorporating a patented hydrogel-based dressing with a breathable outer layer, the dressing can remain adherent even on a patient with hairy or sweaty skin.

The Russell Chest Seal® is a safe and easily applied dressing that can be incorporated into the most basic of first aid measures, with the potential for achieving pronounced clinical improvement of the casualty.

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**Intended use:**

The Russell Chest Seal® is designed for the immediate management of an open chest wound and to minimize the risk of developing tension pneumothorax.

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**Clinical benefits:**

- Valve mechanism allows air to escape from the pleural cavity without being drawn back in.
- Minimizes risk of developing tension pneumothorax and lung collapse.

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

The Russell Chest Seal® is used to cover penetrating chest wounds, and the valve mechanism allows air to escape from the pleural cavity without being drawn back in. This aims to minimize lung collapse associated with this injury and to facilitate lung re-inflation. The device will remain adhered over many hours and is comfortable for the patient in prolonged use.

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