LNCS® TF-I®

SpO2 Reusable Transflectance Sensor








Images	2
en English	3-6
fr French	7-10
de German	11-14
it Italian	15-18
es Spanish	19-22
sv Swedish	23-26
nl Dutch	27-30
da Danish	31-34
pt Portuguese	35-38
zh Chinese	39-42
ja Japanese	43-46
fi Finnish	47-50
no Norwegian	51-54
cs Czech	55-58
hu Hungarian	59-62
pl Polish	63-66
ro Romanian	67-70
sk Slovak	71-74
tr Turkish	75-78
el Greek	79-82
ru Russian	83-86
ko Korean	87-90
ar Arabic	94-91
Performance Specifications	95

LNCS[®] TF-I[®]

SpO₂ Reusable Transflectance Sensor

SpO₂ Reusable Transflectance Sensor

DIRECTIONS FOR USE

Reusable

Not made with natural rubber latex

Non-sterile

Prior to using this Sensor the user should read and understand the Operator's Manual for the Device or Monitor and this Directions for Use.

INDICATIONS

The LNCS® TF-I® Reusable Transflectance sensor is indicated for the continuous non-invasive monitoring of arterial oxygen saturation (SpO₂) and pulse rate for adult and pediatric patients weighing > 30 kg.

CONTRAINDICATIONS

The LNCS TF-I is contraindicated for certain patient positions that affect the monitoring site – see the Warnings Section below. The LNCS TF-I must be removed and repositioned to a different monitoring site at least every two (2) hours. If extended monitoring is required, use of an LNCS single patient adhesive digit sensor is recommended.

DESCRIPTION

The LNCS TF-I sensor is applied to the sensor site using a headband. The LNCS TF-I is for use only with devices containing Masimo SET® oximetry or licensed to use TF-I sensors. Consult individual oximetry system manufacturers for compatibility of particular devices and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor model. The LNCS TF-I series has been verified using Masimo SET Oximetry Technology.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor, cable and sensor before use, otherwise degraded performance and/or patient injury can result.
- The sensor should be free of visible defects, discoloration and damage. If the sensor is discolored or damaged, discontinue use immediately. Never use a damaged sensor or one with exposed electrical circuitry.
- The site must be observed and the sensor repositioned at least every two (2) hours to ensure proper skin integrity and proper alignment.
- Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every (1) hour with poorly perfused patients and move the sensor if there are signs of tissue ischemia.
- During low perfusion, the sensor site needs to be assessed frequently for signs of tissue ischemia, which can lead to pressure necrosis.
- SpO₂ readings may be affected under very low perfusion conditions at the monitored site.
- Do not use tape to secure the sensor to the site; this can restrict blood flow and cause inaccurate readings. Use of additional tape can cause skin damage, and/or pressure necrosis or damage the sensor.
- Misapplied sensors or sensors that become partially dislodged may cause incorrect measurements.
- Sensors applied too tightly may cause erroneously low readings. Applying the headband too tightly can lead to inaccurate saturation measurements, or possibly lead to temporary pressure marks from the sensor.
- Carefully route cable and patient cable to reduce the possibility of patient entanglement or strangulation.
- Do not use this sensor under conditions where the patient is in the Trendelenburg position (head lower than the heart).
- Do not use the sensor during MRI scanning or in a MRI environment.
- Inaccurate readings may be caused by EMI radiation interference.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If sensor is exposed to the radiation, the reading might be inaccurate or not provided for the duration of the active radiation.
- Venous pulsations may cause inaccurate readings (e.g. tricuspid valve regurgitation, Trendelenburg position). Mechanical ventilation can cause increased intra-thoracic pressure and increase venous pooling in the head. This may cause inaccurate readings.
- The pulsations from intra-aortic balloon support can be additive to the pulse rate on the oximeter pulse rate display. Verify patient's pulse rate against the ECG heart rate.
- Venous congestion may cause under reading of actual arterial oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. Trendelenburg position).
- Elevated levels of Methemoglobin (MetHb) will lead to inaccurate SpO₂ measurements.
- Elevated levels of Carboxyhemoglobin (COHb) may lead to inaccurate SpO₂ measurements.

- High levels of COHb or Methb may occur with a seemingly normal SpO₂. When elevated levels of COHb or Methb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
- Elevated Total Bilirubin levels may lead to inaccurate SpO₂ measurements.
- Intravascular dyes such as indocyanine green or methylene blue or externally applied coloring may lead to inaccurate readings.
- Inaccurate SpO₂ readings may be caused by severe anemia, low arterial perfusion or motion artifact.
- Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc. may cause inaccurate SpO₂ readings.
- Inaccurate SpO₂ readings may be caused by elevated levels of dyshemoglobin, hypocapnic or hypercapnic conditions and severe vasoconstriction or hypothermia.
- To prevent interference from ambient light, ensure that the sensor is properly applied. Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.
- High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of the sensor.
- Readings provided with a low signal confidence indicator may not be accurate.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- Do not attempt to remanufacture, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Clean the sensors prior to reuse on multiple patients.
- To prevent damage, do not soak or immerse the connector in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave or ethylene oxide.
- **Caution:** Replace the sensor when a replace sensor message is displayed, or when a low SIQ message is consistently displayed after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- **Note:** The sensor is provided with X-Cal® technology to minimize the risk of inaccurate readings and unanticipated loss of patient monitoring. The sensor will provide up to 8,760 hours of patient monitoring time. Replace the sensor when the patient monitoring time is exhausted.

INSTRUCTIONS

A) Site Selection

- The preferred measuring site is the forehead just above either eye. Do not place the LNCS TF-I on sites with a palpable pulse.
- Site should be cleaned of debris and dry prior to sensor placement.

CAUTION: Before using the sensor, ensure the sensor is physically intact, with no broken or frayed wires or damaged parts.

B) Attaching the sensor to the patient

1. Ensure sensor surface contacting patient skin is dry, clean, free of debris and oil.
2. **OPTIONAL** Refer to **Fig. 1**. For added sensor stability, place one adhesive Pad on sensor surface to contact patient skin, being careful to align the pad such that it does not cover sensor lenses. Once pad is securely placed, remove top release liner, by pulling from release liner tab.

WARNING: Do not use adhesive pads if the patient exhibits allergic reactions to adhesive tape.

3. Refer to **Fig. 2**. Orient the sensor with the emitter and detector 'bumps' toward the patient. Gently press on sensor so that adhesive Pad forms a good contact with patient skin. Orient the sensor so the cable runs upwards and away from the patient.

C) Applying the Headband

CAUTION: Headband must be used to avoid inaccurate readings.

1. Refer to **Fig. 3**. Wrap the Masimo Headband around the patient's head loosely enough to avoid restricting circulation around the site and tightly enough to maintain proper alignment of the sensor.
2. Refer to **Fig. 4**. Secure the Headband using the "hook and loop" tab.
3. Refer to **Fig. 5**. Verify that the sensor is covered and secure behind the Masimo Headband.
4. Use the clothing clip to secure the cable to the patient's gown (not shown).
5. Carefully route the patient cable to reduce the possibility of patient entanglement or strangulation.

D) Attaching the Sensor to the Patient Cable

Refer to **Fig. 6**. Insert the sensor connector completely into the patient cable connector (1). Completely close the protective cover.

E) Disconnecting the Sensor from the Patient Cable

Refer to **Fig. 7**. Lift the protective cover to gain access to the sensor connector (1). Pull firmly on the sensor connector to remove from the patient cable (2). To avoid damage, pull on the sensor connector, not the cable.

CLEANING

CAUTION:

- Do not use undiluted bleach (5–5.25% sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor may occur.
- To prevent damage, do not soak or immerse the connector in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, Gluteraldehyde (Cidex) or ethylene oxide.

To surface clean the sensor

1. Remove the sensor from the patient and disconnect it from the patient cable.
2. Remove Adhesive pad by pulling on white tab and discard. Any residue left by adhesive can easily be removed and cleaned by wiping with 70% isopropyl alcohol.
3. Clean the LNCS TF-I sensor by wiping with a 70% isopropyl alcohol pad. If low-level disinfection is required, wipe with a cloth or gauze pad saturated with a 1:10 bleach / water solution.
4. Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the sensor and cable.
5. Dry the sensor and cable with a clean cloth or dry gauze pad.

To clean or disinfect the sensor using a soaking method

1. Place the sensor in the cleaning solution (1:10 bleach / water solution), such that the sensor pads and desired length of cable are completely immersed.

WARNING: To prevent damage, do not soak or immerse the connector in any liquid solution.

2. Dislodge air bubbles by gently shaking the sensor and cable.
3. Soak the sensor and the cable for at least 10 minutes and not greater than 24 hours.
4. Remove from cleaning solution.
5. Place the sensor and the cable in room temperature sterile or distilled water for 10 minutes. Do not immerse the connector.
6. Remove from the water.
7. Dry the sensor and cable by wiping all surfaces with a clean cloth or dry gauze pad.

SPECIFICATIONS

When used with Masimo SET pulse oximetry monitors, or with licensed Masimo SET pulse oximetry modules and patient cables, the LNCS TF-I sensors have the following specifications:

Sensor:	LNCS TF-I
 Body Weight	> 30 kg
Application Site	Forehead
SpO ₂ Accuracy, No Motion ¹	2 %
Pulse Rate Accuracy, No Motion ²	3 bpm

NOTE: A_{RMS} accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within $\pm A_{RMS}$ of the reference measurements in a controlled study.

¹ The Masimo SET Technology has been validated for no motion SpO₂ accuracy in human blood studies on healthy adult male and female volunteers with light to dark pigmented skin in induced hypoxia studies in the range of 70–100% SpO₂ against a laboratory co-oximeter.

² The Masimo SET Technology has been validated for pulse rate accuracy in bench top testing against a Biotek Index 2 simulator over a range of 25–240 bpm.

COMPATIBILITY

 This sensor is intended for use only with devices containing Masimo SET oximetry or pulse oximetry monitors licensed to use LNCS sensors. Each sensor is designed to operate correctly only on the pulse oximetry systems from the original device manufacturer. Use of this sensor with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of six (6) months. Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ORAL, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER (UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM) EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCT(S) INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

PURCHASE OR POSSESSION OF THIS SENSOR CONFERS NO EXPRESS OR IMPLIED LICENSE TO USE THE SENSOR WITH ANY DEVICE WHICH IS NOT SEPARATELY AUTHORIZED TO USE LNCS SENSORS.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).	Rx ONLY	Federal law (USA) restricts this device to sale by or on the order of a physician
	Consult instructions for use		Lot code		Mark of conformity to European Medical Device Directive 93/42/EEC
	Manufacturer		Catalogue number (model number)		Authorized representative in the European community
	Date of manufacture YYYY-MM-DD		Masimo reference number		Body weight
	Use by YYYY-MM-DD		Storage humidity limitation		Storage temperature range
	Do not re-use/Single patient use only		Do not use if package is damaged		Keep dry
	Pediatric patient		Greater than		Less than
	Non-Sterile		Not made with natural rubber latex		Atmospheric pressure limitation
	Instructions/Directions for Use/Manuals are available in electronic format @ http://www.Masimo.com/TechDocs Note: eIFU is not available for CE mark countries.				

Patents: <http://www.masimo.com/patents.htm>

Masimo, SET, , TF-I, X-Cal and LNCS are federally registered trademarks of Masimo Corporation.

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© 2018 Masimo Corporation
EU Authorized Representative for Masimo Corporation:
MDSS GmbH