next**RESCUE**

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Guangzhou Athmedic Healthcare Co., Ltd.
Room 415, Building 1, No. 9, Shinan second Road,
Dashi Street, Panyu District, 511400 Guangzhou City
P.R. China

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Ohne Naturkautschuklatex hergestellt
Not made with natural rubber latex

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en **INSTRUCTION FOR USE****Intended use:**

Single-use resuscitation face shield to reduce the risk of cross-contamination during mouth-to-mouth resuscitation in adults. Intended for persons trained in the use of resuscitation aids.

Application:

1. Remove the face shield from the pouch. Remove any foreign objects from the patient's airway, if present.
2. Tilt back the head (1). Place the face shield on the patient's face according to the markings ("operator side" readable for the rescuer) (2). Position the valve centrally over the open mouth of the patient. Attach the loops behind the ears.
3. Pinch the patient's nose closed.
4. Breathe slowly through the mouthpiece – check if the chest rises (3).
5. Allow the patient's exhaled air to escape (4).
6. Continue resuscitation until spontaneous breathing resumes.

**Warnings:****Attention!**

1. In case of simultaneous cardiac arrest: Follow current guidelines for cardiopulmonary resuscitation.
2. Does not provide absolute protection but reduces the risk of cross-contamination.
3. Do not use if the packaging is damaged or opened.
4. Not suitable for infants or small children.
5. Do not reuse.

Any serious incident that occurs in relation to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Technical specifications:

- **Material:** Polyethylene (PE)
- **Dimensions (product):** 28.5 cm x 20 cm
- **Dimensions (product packaging):** 5.2 x 5 cm
- **Shelf life:** 5 years from date of manufacture
- **Storage conditions:** -30–50 °C,
Store in a dry, well-ventilated place.
- Protect from sunlight.

Disposal:

Dispose of the product in accordance with the healthcare facility's policies and with local environmental and waste regulations.