REF	Article
20-50	TuBo®-Grip 1,8 x 33 cm + 2,5 x 32,5 cm
20-51	TuBo®-Grip XL 1,8 x 33 cm + 2,5 x 47 cm
20-55	TuBo®-Grip Mittelteil 1,8 x 33 cm
20-52	TuBo®-Grip Nackenpolster 2,5 x 32,5 cm
20-53	TuBo®-Grip Nackenpolster XL 2,5 x 47 cm

Novo Klinik-Service GmbH
Technologiepark West – Haus 15
Zum Frenser Feld 1 | 50127 Bergheim | Germany
Tel. +49 2271 99404-0 | Fax +49 2271 99404-100
info@novo-medde | novo-med.de

IFU_Rev. 03_202408

NOVO®

Gebruiksaanweisung








Instructions for use

TuBo®-Grip

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Endotracheal Tube Holder

TuBo®-Grip / TuBo®-Grip XL

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EN | Instructions for use

Two-part endotracheal tube holder for oralendotracheal tubes, double lumen tubes, laryngealtubes and masks in adolescents and adults.

- Useful life
Change in case of contamination / defect / according to hygienic guidelines. Maximum operating life 24 hours.
- Notice on reporting requirements
Please inform the manufacturer and the competent authority of the Member State in the event of any serious incident that has occurred in relation to the device.

- Notice on medical devices
Medical device of Class I. Accessory product for endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks.

- User training course is recommended before using the product.

- Contraindications
Do not use in case of injuries or accesses in the operating area of the tube holder or allergies to the materials used.

- Warning notices
 - Apply the tube holder in twos and constantly ensure a secure fit of the tube or mask.
 - The safe hold of the fixation has to be checked regularly.
 - If the fixation is too loose, extubation can be caused. A slight tension on the tube holder is optimal.
 - Do not use oils, creams or benzine. Moisture on the tube can negatively affect the hold.
 - The product is a disposable product. The reuse of disposables can result in infection and/or effects on the functioning.

Application by healthcare professionals or users trained by healthcare professionals. Decision on use is made by the treating doctor.

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