READY HEAT THE READY-HEAT FIRST RESPONDER LINE The Temperature Management Solution

The First Responder solution for active patient warming in any casualty situation:



Emergency Rescue

Ready-Heat II

• The #1 selling active warming blanket

blanket, with 8+ hours

• Quick warming

of active warmth

- Ideal for all patient transport and emergency rescue
- Water resistant top layer made of coated medical grade non-woven material and a foam core



Disaster Preparedness

Ready-Heat II Torso

- Compact size and smaller footprint for greater portability
- Deploys in seconds and heats in minutes
- Supplies consistent heat for 8 hours
- 5 year shelf life
- Water resistant top layer made of coated medical grade non-woven material and a foam core



Air Transport

Ready-Heat Infant Cocoon

• Designed to actively warm the most precious cargo

- Consists of 6 thermal pads strategically placed to warm infant's body
- Hood assists in the retention of body warmth and the padded foam core and dual hook and loop design provide maximum comfort and fit









For the First Responder

Stay warm on the job with proprietary Ready-Heat technology in a disposable heated vest.

Include Ready-Heat in Your:

- First Aid Kits
- Wilderness Survival Kits
- Mountain Rescue Kits
- Pediatric Transport Kits
- Water Rescue Kits
- Trauma Kits





Visit our website for more detailed information and IFUs - www.ready-heat.com/IFU www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130

0459 ISO 13485:2016

RHB-3 Rev1 July 2021

wa
car
• Co
рас
pla

READY PRODUCT HEAT CHART

4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8 Warming Temperature:

READY-HEAT II BLANKET

104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 48" / 86 cm x 122 cm

Weight: 1.75 lbs. / 794 gm

Part Number / Per Case: G12RH2 / Case of 8

Warming Temperature:

100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours

Size: 34" x 48" / 86 cm x 122 cm

Description: First Responder Blanket. Used

to Treat Shock, Trauma and Pre-Hypothermia.

Has a foam core for water resistance and is

lightweight and comfortable. Quick to reach

operating temperature. Works well at high

Weight: 1.5 lbs. / 794 gm

Description: Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

NSN #: 6532-01-538-0759

6 PANEL BLAI	NKET NSN #: 6532-01-525-4063	
Part Number / Per Case: S6RHLG / Case of 8		
	Warming Temperature: 104°F / 40°C in 15-20 Minutes	
	Maintains Heat: Up to 10 Hours	
	Size: 34" x 60" / 86 cm x 152 cm	
	Weight: 2.5 lbs. / 1.13 kg	

Description: Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

READY-HEAT II TORSO BLANKET

Part Number / Per Case: G6RH2T / Case of 12

Warming Temperature: 100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours

Size: 34" x 28" / 86 cm x 61 cm

Weight: 14 oz. / 397 gm

Description: First Responder Blanket. Used to treat shock, trauma and pre-hypothermia. Smaller overall footprint both in bag and on patient. All the benefits of the full-sized blanket. G12RH2.

TEMPERATURE MANAGEMENT FULL BODY

ſ

Part Number / Per Case: SB9RH9120 / Case of 10		
	Warming Temperature: 104°F / 40°C in 15-20 Minutes	
	Maintains Heat: Up to 10 Hours	
	Size: 34" x 60" / 86 cm x 152 cm	
	Weight: 18 oz. / 510 gm	
	Description: Patient Warming Blanket. Full mobile, can move with the patient from pre-op to intra-op and post-op. Prevents	

1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



104°F / 40°C in 15-20 Minutes

Warming Temperature:

Maintains Heat: Up to 10 Hours

Size: 16" x 18" / 41 cm x 46 cm

Weight: 7 oz. / 198 gm

Description: Mini Blanket. Suitable for injured extremities. Fully adustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

READY-HEAT PERSONAL USE VEST

Part Number / Per Case: GRHV06A / Case of 12



Warming Temperature: 100°F / 38°C in 10-15 Minutes

Maintains Heat: Up to 8 Hours

Size: Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

Weight: 10 oz. / 284 gm

Description: Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

INFANT WARMING COCOON

altitudes.

Part Number / Per Case: GIW6C / Case of 12



Warming Temperature: 100°F / 38°C in 15 Minutes

Maintains Heat: Up to 8 Hours

Size: 31.5" x 28" / 80 cm x 71 cm

Weight: 12 oz. / 340 gm

Description: Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

lν hypothermia, increases patient comfort for higher quality of care, effective solution for infection control.

TEMPERATURE MANAGEMENT HALF BODY

Part Number / Per Case: SB6RH9120 / Case of 20

Warming Temperature: 104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 24" / 86 cm x 61 cm

Weight: 12 oz. / 350 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from preop to intra-op and post-op. Half-sized blanket can be moved and adjusted for different types of surgeries or access. All the benefits of the full-sized blanket, SB9RH9120.

MANUFACTURED BY TECHTRADE LLC IN THE USA

Visit our website for more detailed information and IFUs - www.readv-heat.com/IFU www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130