

SX RED KEY

Pick-up stretcher



SX RED KEY pick-up stretchers are lifting devices designed to move severely injured patients to other transport devices.

Specific features

- Extremely light and resistant structure
- The possibility to adjust the length confers high versatility
- Foldable to save space when stored
- Smooth surfaces easy to clean
- The numerous handles make the transport easier
- Belts with transport bag

Standard equipment

- K Belt 3 – Three one piece belts STX 597

Accessories

- | | |
|----------|-------------------------------------|
| ST05009A | SX Bag – Orange transport bag |
| ST05018C | SXV – 20 G Vertical fixing system |
| ST05019C | SXO – 20 G Horizontal fixing system |

Compliant with

UNI EN 1865-1

Technical data ¹

Shortened length	1660 mm
Intermediate length 1	1780 mm
Intermediate length 2	1900 mm
Maximum length	2050 mm
Length when folded	1190 mm
Thickness when folded	70 mm
Thickness when extended	60 mm
Width	440 mm
Materials	Aluminium
Load capacity	170 kg
Weight	7,6 ± 0,4 kg

¹ Length measures are subjected to tolerances of ± 10 mm

Class I MD compliant with UE Reg. 2017/745

SPENCER SXR PICK UP STRETCHER,
2 BLADES W/SET 3 BELTS

ST05002A

CND Classification V9099

Registration number 91818



Rev.0 (07/10/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it