

Product Data Sheet

Polycarbonate Eye Shield (PES)

Description: Convex, rigid eye protection made of transparent, orange polycarbonate. Shatterproof and shatterproof. With ventilation holes for air circulation. For covering and securing eye trauma. Complies with the care recommendation of the TCCC guidelines. Protects the injured eye and prevents external pressure being exerted on the injury. Can be used for both eyes.

REF: NR85-ES-PCO



Product features	
Material	Polycarbonate
Color	Orange
Dimensions (basic package)	7,6 x 6,2 cm (tolerance ±0,3cm)
Weight (incl. packaging)	Approx. 7,5 g

Regulatory and market-specific information	
Regulation	CE-compliant medical device class is according to Regulation (EU) 2017/745 (MDR)
NSN	6515-01-647-5456
Country of manufacture	CN
Customs tariff number	39269097900

Packaging properties, storage	
Basic package	Content: 1 piece PE pressure seal bag (outside) with paper insert
Sterile product	No
Shelf life	4 years from date of manufacture
Collective pack	Box with 100 pieces
Shipping packaging	Carton of 10 collective packs = 1000 pieces
Shipping packaging:: Dimensions, weight	Approx. 47 x 29 x 33 cm, Approx. 9,4 kg
Storage information	Protect from sunlight Store in a dry place

European Importer:
nextRESCUE GmbH
Brennhorster Str. 12
32479 Hille, Germany

EC REP Mednet EC-REP GmbH
Borkstrasse 10
D-48163 Münster

Anping Guardian Medical Equipment
Co., Ltd. 500 Meters, Nan Sucun, Anping
County Hebei, China 053600