

STX 540

Fixed lifting bridle system



These harnesses are accessories for basket stretchers, to be used solely for lifting from a floor-mounted position.

Specific features

- 4 straps to be fixed by means of a noose knot
- Reinforced seams
- Carabiner compliant with EN 362 and EN 365

Technical data

Belt length	1600/1550 mm
Belt width	30 mm
Weight	0.8 ± 0.1 kg
Material	Nylon, aluminium, steel
Fastening type	Loop
Number of fastening bands	4
Compatibility	All basket stretchers manufactured by Spencer
Maximum applicable static safety load	450 kg (1000 lbs)

Class I MD compliant with UE Reg. 2017/745

STX 540 FIXED LIFTING
BRIDL.-LIFTING BRIDLE SYST.

ST04522B

CND Classification



V0899

Rev.0 (27/04/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it