

FCP-09, 8" Emergency Bandage®

aka The Israeli Bandage®

- Best-in-class multiple use bandage addressing eviscerations, large area wounds, and amputations
- 12"x12" non-adherent pad attached to an 8" wide elastic leader
- Direct pressure can either be used or bypassed depending on wound type
- 8-year shelf life



Extended Coverage for Major Bleeding

The FCP-09 Emergency Bandage® provides extensive coverage for eviscerations and the ability to provide pressure for large area wounds and amputations.

The versatile design is an effective solution for the control of severe bleeding, eviscerations, large area wound coverage, splinting, and more. The device is intuitive and simple to use for trained and untrained responders.



8" Green Emergency Bandage®

Abdominal/Large Area Wounds



Product Highlights

Large abdominal Emergency Bandage® design

Intuitive design facilitates easy, rapid application


Multiple use bandage with extensive wound coverage

This extended coverage all-in-one bandage is an ideal component for any medical kit

- Bleeding control, compression, and coverage for large area wounds and amputations
- Extended coverage and containment for eviscerations
- Easy, intuitive application
- Integrated pressure bar may be used to apply direct pressure to amputations and large area wounds
- Integrated pressure bar may be bypassed to allow wound coverage without applying pressure on eviscerations
- Large single sterile pad
- Elasticized wrapping leader
- Stop-N-Go™ release for controlled wrapping
- Simple closure bar

Specifications

Description	8" (20cm) wide abdominal/large wound bandage
Part Number	FCP-09
NSN	6510-01-532-6656
Pad Size	12" wide x 12" long (30cm x 30cm)
Full Length	75" (190cm) unstretched, including pad
Shelf Life	8 years
Colors	OD green elastic leader, white pad



PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

