




BloodCOMM

REMOTE BLOOD PROGRAM MANAGEMENT SYSTEM

BloodCOMM is a highly intuitive blood program management software application that allows you to wirelessly connect to existing APRU 6L Portable Medical Blood Refrigerators and Delta ICE 2L Smart Blood Coolers. It is the next evolution in mobile blood management that provides accountability and ensures compliance with Blood Banks, Hospitals and Agency Directors.

BloodCOMM is accessible on the app or through your web browser at <https://blood.dev/>.



- Real-time temperature monitoring and recording
- Close-range (Bluetooth) and Remote (Cellular) Monitoring
- View temperature history and detailed graphs
- FREE Bluetooth capability for Delta ICE and APRU
- Alerts of temperature excursions, approaching excursions, low battery
- Download, export and email reports in PDF or CSV
- Unlimited users per organization
- Assign custom nicknames per device
- Affordable cellular subscription plans
- Paid subscription required for remote, cellular connection for Delta ICE & APRU



PRO

*annual fee/device applies

Remote Accountability

- Includes:
- Bluetooth Features
 - Cellular Connectivity
 - Remote Monitoring
 - Remote Reports
 - Remote Support
 - Push Notifications
 - Unlimited Team Size



PREMIUM

*annual fee/device applies

Advanced Features

- Includes:
- PRO Features
 - 24/7 Technical Support
 - Advanced Program Metrics
 - Blood Bag Scanning
 - Metrics Dashboard
 - Cold Chain Annotations

**Not yet available*

SUBSCRIPTION PLANS



DELTA
DEVELOPMENT TEAM

sales@deltadevteam.com

deltadevteam.com

[@deltadevelopmentteam](https://t.me/deltadevelopmentteam)