

Combat CARL Care and Evacuation Manikin

Item no.: T3C




*Experts in medical education
since 1950*



ERLER ZIMMER

Erler-Zimmer GmbH & Co. KG · Hauptstr. 27 · D-77886 Lauf
Tel.: (+49) 07841/67191-0 · info@erler-zimmer.de · www.erler-zimmer.de



Combat CARL offers all functions for TCCC and TECC training.

Life-size (183cm), 75 kg, with highly mobile extremities, no sensitive electronics, no gimmicks: Ideal for treatment and transport in extreme conditions. Suitable for practicing the standard algorithms.

MARCH

- M** Three included deep wounds can be applied to 7 sites and connected to bleeding modules. Wound packing is possible for all wounds, as is the application of tourniquets.
- A** Airway management with oro- and nasopharyngeal as well as supraglottic devices.
- R** Thoracic relief in monaldi and Bülau position
- C** Hemostasis and access to the vascular system: i.v. in the forearm, i.o. sternal, humerus head and proximal tibia with multiple-use pads.
- H** Voice output from prerecorded files or directly.

Robust and insensitive to water, dirt, cold and heat.

- Thorax relief Monaldi
- Thorax relief Bülau
- IV access possibility

- I.O Access FAST
- I.O Access
- Exchangeable eyes with normal, constricted or dilated pupil



- Voice output via loudspeaker in the head by means of sound files or via microphone function
- Airway management oro - and nasopharyngeal

- Openings for wounds, IM puncture pads and blind covers
- Robust and insensitive to water, dirt and other adverse environmental conditions



**Airway management
nasopharyngeal**



**Airway management
supraglottisch**



Patient voice



**Exchangeable
eyes**



**Wounds in the
thoracic region**



**Relief puncture according
to Monaldi**




**Sternal intraosseous
approach**



Wound-Packing




i.V. additions



Chest Drainage



**IO access
Humerus**



**IO access
Tibia**



TQ - placement



Bandaging techniques



Erler-Zimmer GmbH & Co.KG

Hauptstraße 27

D-77886 Lauf

Tel.: (+49) 07841/67191-0

Fax: (+49) 07841/67191-99

E-Mail: info@erler-zimmer.de

Available 24/7:
www.erler-zimmer.de

Follow us!

