



CERTIFIED
MEDICAL
DEVICE



STAIR
CLIMBING AND
DESCENDING



LOADING
CAPACITY



BATTERY
AUTONOMY



POWERED
MECHANISM

GLIDE UP AND DOWN
STAIRS WITH THE PUSH
OF A BUTTON



A FULLY POWERED
SYSTEM FOR SECURE
STAIR ASCENTS AND
DESCENTS

NEXTEP

Advanced Powered Chair for Effortless Stair
Ascending and Descending
Ensuring Safe Transport for the Mobility
Impaired



 SPENCER

NEXTEP

REVOLUTIONIZING ACCESSIBILITY WITH OUR
MOTORIZED STAIR CLIMBING CHAIR.

Designed to empower individuals with limited mobility, this medical device ensures safe, controlled, and effortless navigation of stairs, both up and down.

Equipped with innovative powered tracks, it offers unparalleled stability and ease of use. Enhanced with four swivel castors and an efficient quick-closing system, NEXTEP boasts exceptional maneuverability, making it an indispensable tool for enhancing mobility and independence.



- ✓ EN ISO 9001:2015
- ✓ EN 1789
- ✓ 2 YEARS WARRANTY

ELEVATE MOBILITY WITH EFFORTLESS STAIR NAVIGATION.

Intuitive display with battery charge status indication and controls to vary motor speed and activate LED lights.

The telescopic handle is adjustable in height, making it easy to use.

Convenient and easy to use, NEXTEP requires no effort when descending or ascending stairs.

Made of aluminium alloy with black ABS seat and backrest.

Rubber tracks provide excellent grip on most surfaces.

2 large diameter rear wheels with brakes facilitate transportation on flat surfaces and uneven ground.

The front wheels swivel for easy steering and manoeuvrability of the chair in tight spaces.



**SAFELY EMPOWER MOBILITY
FOR ALL: USER-FRIENDLY
SOLUTIONS FOR EVERY
CAREGIVER.**

Easy to fold and stow in the
car trunk or in a van.



NEXTEP STORAGE SPACE
115 X 50 X 27 cm



Recommended use by
2 operators

MINIMUM LENGTH	80 cm
MAXIMUM LENGTH	110 cm
WIDTH	50 cm
MINIMUM HEIGHT	110 cm
MAXIMUM HEIGHT	160 cm
WEIGHT	30 kg
LOAD CAPACITY	180 kg
MATERIAL	Aluminium alloy, ABS
MOTOR	24V-200W or 36V-300W
BATTERY LIFE	3-5 years
SPEED	45 steps per min.
FULL CHARGING TIME	5-6 hours



**SETTING THE
STANDARD
IN MEDICAL
TRANSPORT:
YOUR GLOBAL
PARTNER IN
EMERGENCY
INNOVATION**

As global leaders in emergency and medical transport equipment, we proudly stand among the largest manufacturers in the field. Our portfolio encompasses over 1,700 premium Italian-made products, spanning patient recovery and transport, immobilization, oxygen therapy, and first aid.

Established in 1989, we are committed to engineering devices that excel in critical moments—ensuring unparalleled performance and safety when it matters most.

Class I Medical Device



SPENCER ITALIA SRL

Via Provinciale 12, 43038, Sala Baganza (PR) Italy
+39 0521 541111 | export@spencer.it