

Product Data Sheet

Bow Emergency Trauma Bandage

Description: Trauma emergency bandage type „Israeli Bandage“, width 10 cm, highly absorbent wound dressing 10 x 20 cm made of non-adherent material, bandage length approx. 4.5 m (stretched), with integrated pressure applicator in arch form and a closure clip for secure fixation of the bandage. Roll-off protection for controlled application. Sterile and vacuum packed.

REF: NR70-VSIB-4-01



Product features

Material	Blend of cotton, spandex and polyester. Plastic parts made of polycarbonate.
Color	Green
Dimensions (basic package)	Ca. 12 x 8 x 4 cm
Weight (incl. packaging)	Ca. 70 g

Regulatory and market-specific information

Regulation	CE-compliant medical device class Is according to Regulation (EU) 2017/745 (MDR)
NSN	6510-12-376-4655
Country of manufacture	CN
Customs tariff number	30059099

Packaging properties, storage

Basic package	Content: 1 piece Double-layer packaging made of PE film (outside) and nylon film (inside)
Sterile product	Yes
Shelf life	8 years from date of manufacture
Collective pack	/
Shipping packaging	200 pieces
Shipping packaging:: Dimensions, weight	Approx. 55 x 35 x 45 cm, Approx. 12 kg
Storage information	Protect from sunlight Store in a dry place



European Importer:
nextRESCUE GmbH
Brennhorster Str. 12
32479 Hille, Germany



Caretechion GmbH
Niederrheinstr. 71
40474 Düsseldorf
Germany



Anji Hengfeng Sanitary Material Co., Ltd.
Ancheng, Dipu Town, Anji County
313300 Zhejiang, P.R. China