



Product Data Sheet en

Red Point Gaze Compressed, z-fold

Description: White gauze bandage made of 100 % cotton, 11.43 cm x 3.75 m (4.5" x 4.1 yds), 6-ply, compressed and vacuum-packed, sterile, high absorbency, double-layered packaging for easy and hygienic handling, wide range of applications for wound care.

REF: NR70-VSCGF-01



RapidColour System®
Dressings
Tactical Gray

Product features

Material	100 % cotton
Color	White
Dimensions (basic package)	Approx. 15 x 10 x 1,5 cm
Weight (incl. packaging)	Approx. 40 g

Regulatory and market-specific information

Regulation	CE-compliant medical device class Is according to Regulation (EU) 2017/745 (MDR)
Country of manufacture	CN
Customs tariff number	3005901000
UDI	06971519611091
Manufacturer-Reference	HF-CG012

Packaging properties, storage

Basic package	Content: 1 piece Double-layer packaging made of PE film (outside) and nylon film (inside)
Sterile product	Yes
Shelf life	8 years from date of manufacture
Collective pack	/
Shipping packaging	200 pieces
Shipping packaging: Dimensions, weight	Approx. 55 x 35 x 20 cm, Approx. 9 kg
Storage information	Protect from sunlight. Store in a dry place.



Europäischer Importeur:
nextRESCUE GmbH
Brennhorster Str. 12
32479 Hille
Germany

EU REP

Caretechion GmbH
Niederrheinstr. 71
40474 Düsseldorf
Germany

CH REP

SUNGO Technical
Service GmbH
Bahnhofstrasse 21
6300 Zug
Switzerland



Anji Hengfeng Sanitary
Material Co., Ltd. No. 128,
Liuguanli Road, Dipu Sub-
district, Anji County,
313300 Zhejiang, China