


SICH

T O U R N I Q U E T



SICH Ukraine is a leading Ukrainian manufacturer of **SICH Tourniquets** — the **No. 1** tool in modern civilian and military first aid kits.

Since 2016 and to this day, the official supplier to the Ministry of Defense of Ukraine, the National Guard, and the National Police.



The SICh Tourniquet is a **next-generation hemorrhage** control device, developed with insights from modern armed conflicts to save the lives of military personnel, first responders, police officers, drivers, travelers, and civilians — in any climate conditions.



ICE C S



Combat-proven



Functions when wet, bloody, dirty, snow-covered, in freezing cold or extreme heat



Withstands extreme stress and pressure



Can be applied with **one hand**



in under **10 seconds**



Non-reflective surface — no light glare



The **SICH Tourniquet** features a **specialized design** that works reliably in all weather conditions and stops massive limb bleeding within seconds.




UNIQUE ADVANTAGES

THAT OUTPERFORM COMPETITORS




METAL HARDWARE

Lightweight yet durable — made of **aircraft-grade aluminum** to withstand extreme physical stress.




ZERO TRAUMA

Soft coating, a **solid one-piece buckle** and **platform**, and a patented internal **telescopic mechanism** eliminate the risk of pinching soft tissues or causing tissue atrophy during application.



FIXATION TRIANGLE

A unique, **patented triangular design** ensures fast and secure fixation of the windlass and the free strap end.



ERGONOMICS

Not a single unnecessary millimeter.

The platform and windlass are optimally sized — refined through extensive field testing under extreme conditions. Delivers maximum surface pressure to stop blood flow more effectively than other tourniquets.

TIME STRAP FOR MARKING




A special **gray coating** allows clear marking even without a marker — using any sharp object (e.g., a cartridge, nail, or similar).



ULTRA-STRONG BONDING TECHNOLOGY

A **patented bonding technology** ensures exceptionally **strong adhesion of the strap**, guaranteeing reliable performance when wet, dirty, snow-covered, in freezing cold or extreme heat.




The only tourniquet of its kind capable of stopping bleeding on limbs with diameters **from 5 to 40 cm**, without pinching soft tissues, and precisely engineered so every millimeter works efficiently.

EFFECTIVE ON SLENDER LIMBS AND K9 UNITS

Specifically designed to work even on the thinnest human limbs, and can be used on service dogs (K9).





Today, the SICh Tourniquet is the most widely used hemorrhage control device in Ukraine and **has already saved thousands of lives.**

SCH

