

Product: TRIANGULAR BANDAGE

Reference code: 1889



1. General description

Two individually wrapped triangular bandages packed in a printed box together with safety-pins.

2. Intended use

The Triangular Bandage can be used as a sling to relieve pressure on a wounded arm or tied around the head as an outer covering for a dressing or wound pad.

3. Technical description

Product dimensions (l x w x h) cm:

Triangular bandage 96 x 96 x 136
 Package 12,5 x 3,5 x 6,5

Material

Bandage White spunbonded polypropylene

4. Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
51618330	Triangular bandage	2
51618760	Safety-pins	6
51616580	Printed box	1

Transport package

Number of consumer packages in transport package: 24
 Dimensions (l x w x h) cm: 26 x 21 x 18
 Weight (kg): 1,3

5. Labeling

On consumer package

Cederroth address, telephone- and faxnumber, web-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Symbols: LOT, USE BY

On transport package

Cederroth address, telephone- and faxnumber, web-address
 Product name, Reference code
 Lot-number, Expiry date, EAN-code
 Quantity
 Symbols: LOT, USE BY



PRODUCT SPECIFICATION

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years based on production date

Authorized by: Teija Ålander
Quality Assurance & Regulatory Manager

Date: 2009-05-05



www.firstaid.cederroth.com

Cederroth Industrial Products AB, Box 715, Kanalvägen 10A, SE 194 27 Upplands Väsby, Sweden.
Phone: +46 8 590 96 300, Fax: +46 8 590 30 953, E-mail: firstaid@cederroth.com