

TEMPWATCH

Delta Development's Tempwatch is designed to interface with the Delta ICE 2L PCM Pack during conditioning to monitor the embedded temperature sensors and alert the user with a unique alarm melody once the PCM pack has reached 2°C. The safe temperature range for blood storage is between 2-6°C.



BENEFITS

- Set and Forget - Automatically Conditions the PCM Pack, taking out any Guesswork
- Displays the Current Status of PCM Pack (Too Cold, In-Range, Too Warm, Error)
- Automatically Alerts the User once the Pack is Conditioned to 2°C
- Maximizes the Holdover Time for Each Conditioned PCM Pack
- Increases Efficiency

SPECIFICATIONS

- Water Resistant
- Designed for Bench-Top Operation
- Interfaces with ICE 2L PCM Pack
- Dimensions: 10"L x 10"W x 2"H
- Unit Weight: 1 lb
- Operating Temperatures: -20 - 50°C
- Storage Temperatures: -20 - 70°C
- USBC Cable and AC/DC Adapter Included

