Technical Data Sheet

Rescue Trade Bed-Set Basic UDI-DI: 426072773RT-BEDSET7J

Rescue Trade Disposable Blanket – 500 g

Item: RTEDPWB500

Description and Application

The disposable patient blanket is intended for short-term single use in healthcare facilities, such as hospitals, ambulances, patient transport and nursing homes, and home care, and provides a physical barrier against the direct transmission of germs infectious potentially to the patient. It also increases comfort and well-being during the patient's stay or during transport.

Material and Characteristic

- Color: blue
- Filling: Polyester 180g/m²
- Outer Layer: 2 Layers PP-nonwoven, 30g/m²
- Weight: approx. 500 grams
- Size: approx. 190x110 cm
- Individually hygienically and space-savingly packed in polybags

Packaging Information

• Unit:

🐨 = 40 pcs/ctn = 20 ctn/Pallet

- Carton Size: approx. 60x40x40 cm
- Net Weight Ctn: approx. 20 kg
- Gross Weight Ctn: approx. 21 kg

Additional Information

MD Medical Device Class I according to EU Regulation 2017/745 MDR

• GTIN Carton: 4260727730123

- GTIN Pieces: 4260727730680
- GMDN-Number: 60823
- FDA-CFR Code: §880.6060



Rescue Trade OHG Am Bahndamm 7 DE-63683 Ortenberg SRN-Number: DE-MF-000008980 FDA-Registration-Nr.: 3027240400

V5 Date 27.11.2023 No guarantee for misprints or mistakes.

Contact:









