SMRT™ TASK TRAINER - HIP WITH LARGE CALIBER GUNSHOT WOUND

Learn more about this and other *Products with a Mission®* at *www.NARescue.com*

STRATEGIC MEDICAL RESEARCH & TRAINING™

Large caliber rifle wounds with large cavities create significant challenges. The Hip with Large Caliber Gunshot Wound is a unique wound captures the hemorrhage point identification wound packing challenges unique to high velocity wounds.

This large cavity wound simulates a high-velocity gunshot wound to the hip with exit through the lateral thigh with extensive undermined tissue. This wound creates a unique assessment and treatment challenge as the entrance wound is often initially missed and the source of the bleeding must be accessed through the wound cavity. This wound may also be used for blast scenarios due to the large irregular wound pattern.

The innovative and wearable SMRT™ Moulage Task Trainers from Strategic Medical Research & Training™ create an unmatched experience of realism featuring realistic moulage and active bleeding across a wide spectrum of woundings and interventions. Made with durable materials and designed for a "one size fits most" range of limb and body sizes, the SMRT™ Moulage Task Trainers are engineered to work with live role players as well as to augment existing simulation systems—while allowing flexibility in placement and utilization. All SMRT™ Moulage Task Trainers work with SMRT™ Trauma Tables, SMRT™ Bleeding Systems, and common hydration systems, making them standalone and ready to use, while also allowing ease of integration into existing training systems. Utilized by Military, Law Enforcement, EMS, and Hospitals to create unsurpassed realism in hemorrhage control training, SMRT™ Training Aids are realistic, cost effective, highly mobile, durable, and low maintenance.



FEATURES/SPECIFICATIONS

- Unsurpassed realism in hemorrhage control training
- One size fits most range of limb and body sizes
- Hip with Large Caliber Gunshot Wound compatible with other SMRT™ Task Trainers
- Available in Light Tone or Deep Tone

SMRT™ HIP W/ GSW - LARGE CALIBER

ITEM #	COLOR
93-1199	Light Tone
93-1200	Deep Tone

** Explore other SMRT™ moulage task trainers at www.NARescue.com*

DIMENSIONS

- Packaged Dimensions: L 11 in. x W 18 in. x D 4 in.
- Weight: 5.1 lbs.

NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277